

Safety Data Sheet

according to UK REACH Regulation

Identium Medium

Revision date: 13.09.2022

Product code: 001166

Page 1 of 6

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Identium Medium

Further trade names

Identium Medium

Identium Medium Fast

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture

Preparation of silicone polymers and fillers with platinum-catalyst (ppm-range).

1.3. Details of the supplier of the safety data sheet

Company name:	Kettenbach GmbH & Co. KG	
Street:	Im Heerfeld 7	
Place:	D-35713 Eschenburg	
Telephone:	+ 49(0) 2774/705-0	Telefax: +49 (0)2774/705-66
e-mail:	info@kettenbach.com	
Internet:	www.kettenbach.com	
Responsible Department:	Regulatory Affairs	

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

GB CLP Regulation

This mixture is not classified as hazardous in accordance with GB CLP Regulation.

2.2. Label elements

GB CLP Regulation

Precautionary statements

P501	Dispose of contents/container to Refer to any Community provisions relating to waste..
P305+P354+P338	IF IN EYES: Immediately rinse with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P280	Wear protective gloves.

Additional advice on labelling

There is no requirement for the product to be specially labelled according to EC directives or the corresponding national laws.

2.3. Other hazards

No risks worthy of mention.

SECTION 3: Composition/information on ingredients

3.2. Mixtures

Chemical characterization

Preparation of silicone polymers and fillers with platinum-catalyst (ppm-range).

Safety Data Sheet

according to UK REACH Regulation

Identium Medium

Revision date: 13.09.2022

Product code: 001166

Page 2 of 6

Hazardous components

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification (GB CLP Regulation)			
63148-57-2	4-Methoxyphenol			< 5 %
	Repr. 2; H361f			
68478-92-2	Platinum			<2 %
	270-844-4			

Full text of H and EUH statements: see section 16.

SECTION 4: First aid measures

4.1. Description of first aid measures

After inhalation

not applicable

After contact with skin

In general the product is not skin irritative.

After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water. In case of troubles or persistent symptoms, consult an ophthalmologist.

After ingestion

Seek medical attention if problems persist.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

Water fog. Foam. Carbon dioxide (CO₂).

5.2. Special hazards arising from the substance or mixture

Not applicable.

5.3. Advice for firefighters

Co-ordinate fire-fighting measures to the fire surroundings.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General advice

Not required.

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up

Other information

Take up mechanically.

6.4. Reference to other sections

No release of hazardous substances.

SECTION 7: Handling and storage

Safety Data Sheet

according to UK REACH Regulation

Identium Medium

Revision date: 13.09.2022

Product code: 001166

Page 3 of 6

7.1. Precautions for safe handling

Advice on safe handling

No special precautionary measures are necessary.

Advice on protection against fire and explosion

No special precautionary measures are necessary.

Advice on general occupational hygiene

The usual precautions when handling chemicals have to be obtained.
Avoid contact with skin, eyes and clothes.

Further information on handling

No special precautionary measures are necessary.

7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels

No special requirements.

Hints on joint storage

Not required.

Further information on storage conditions

No special requirements.

7.3. Specific end use(s)

Not required.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Exposure limits (EH40)

CAS No	Substance	ppm	mg/m ³	fibres/ml	Category	Origin
56-81-5	Glycerol, mist	-	10		TWA (8 h)	WEL

8.2. Exposure controls

Appropriate engineering controls

not applicable

Individual protection measures, such as personal protective equipment

Eye/face protection

Not required.

Hand protection

Single-use gloves.

Suitable material: PE (polyethylene).

Unsuitable material: NR (Natural rubber (Caoutchouc), Natural latex).

For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

Skin protection

Not required.

Respiratory protection

With correct and proper use, and under normal conditions, breathing protection is not required.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state:

Paste

Safety Data Sheet

according to UK REACH Regulation

Identium Medium

Revision date: 13.09.2022

Product code: 001166

Page 4 of 6

Colour: Catalyst: grey-white. Base: violet.
Odour: odourless

Changes in the physical state

Melting point/freezing point: not applicable

Boiling point or initial boiling point and boiling range: not applicable

Explosive properties

not Explosive.

pH-Value: neutral

Water solubility: insoluble

Density (at 20 °C): ca. 1,5 g/cm³

SECTION 10: Stability and reactivity

10.1. Reactivity

None by intended use.

10.2. Chemical stability

None by intended use.

10.3. Possibility of hazardous reactions

None by intended use.

10.4. Conditions to avoid

None by intended use.

10.5. Incompatible materials

None by intended use.

10.6. Hazardous decomposition products

None by intended use.

Further information

Please pay attention to the Directions for Use.

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in GB CLP Regulation

Acute toxicity

Acute toxicity, oral LD50: > 2000 mg/kg Rat

Irritation and corrosivity

Irritant effect on the skin: non-irritant.

Irritant effect on the eye: non-irritant.

Sensitising effects

no danger of sensitization.

11.2. Information on other hazards

Other information

None by intended use.

SECTION 12: Ecological information

12.1. Toxicity

product has not been tested.

12.2. Persistence and degradability

product has not been tested.

12.3. Bioaccumulative potential

Safety Data Sheet

according to UK REACH Regulation

Identium Medium

Revision date: 13.09.2022

Product code: 001166

Page 5 of 6

product has not been tested.

12.4. Mobility in soil

product has not been tested.

12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to UK REACH.

12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

Further information

Do not allow to enter into surface water or drains.

1 (slightly hazardous to water)

Doesn't get into the sewage water as long as the process is carried out according to regulations.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal recommendations

Remove according to the regulations.

List of Wastes Code - residues/unused products

070217 WASTES FROM ORGANIC CHEMICAL PROCESSES; wastes from the MFSU of plastics, synthetic rubber and man-made fibres; waste containing silicones other than those mentioned in 07 02 16

List of Wastes Code - used product

200399 MUNICIPAL WASTES (HOUSEHOLD WASTE AND SIMILAR COMMERCIAL, INDUSTRIAL AND INSTITUTIONAL WASTES) INCLUDING SEPARATELY COLLECTED FRACTIONS; other municipal wastes; municipal wastes not otherwise specified

Contaminated packaging

Completely emptied packings can be re-cycled.

SECTION 14: Transport information

Land transport (ADR/RID)

14.1. UN number or ID number: N/A

14.2. UN proper shipping name: N/A

14.3. Transport hazard class(es): N/A

14.4. Packing group: N/A

Inland waterways transport (ADN)

14.1. UN number or ID number: N/A

14.2. UN proper shipping name: N/A

14.3. Transport hazard class(es): N/A

14.4. Packing group: N/A

Marine transport (IMDG)

14.1. UN number or ID number: N/A

14.2. UN proper shipping name: N/A

14.3. Transport hazard class(es): N/A

14.4. Packing group: N/A

Air transport (ICAO-TI/IATA-DGR)

14.1. UN number or ID number: N/A

14.2. UN proper shipping name: N/A

14.3. Transport hazard class(es): N/A

14.4. Packing group: N/A

Danger releasing substance: N/A

Safety Data Sheet

according to UK REACH Regulation

Identium Medium

Revision date: 13.09.2022

Product code: 001166

Page 6 of 6

14.6. Special precautions for user

N/A

14.7. Maritime transport in bulk according to IMO instruments

N/A

Other applicable information

No dangerous good in sense of these transport regulations.

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

National regulatory information

Water hazard class (D): 1 - slightly hazardous to water

15.2. Chemical safety assessment

For this substance a chemical safety assessment is not required.

SECTION 16: Other information

Changes

Update 2022 09 13

Relevant H and EUH statements (number and full text)

H361f Suspected of damaging fertility.

Further Information

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights.

The receiver of our product is singularly responsible for adhering to existing laws and regulations.

(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)