



Safety Data Sheet Permlastic Base (Regular Body, Light Body & Heavy Body)

SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

1.1. Product identifier

Product name : Permlastic Base (Regular Body, Light Body & Heavy Body)

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses

Main use category : Professional use

Function or use category : Dental materials

Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Kerr Italia S.r.l.
Via Passanti, 332
84018 Scafati (SA) - Italy
T +39-081-850-8311
E-mail: safety@kerrhawe.com

Manufacturer

Kerr Italia S.r.l.
Via Passanti, 332
84018 Scafati (SA) - Italy
T +39-081-850-8311
E-mail: safety@kerrhawe.com

Contact person : safety@kerrhawe.com - tel. 00-800-41-050-505 (08.00-17.00)

1.4. Emergency telephone number

Emergency number : CHEMTREC® Emergency Call Center. Emergency Telephone Number (for USA only) 001-800-424-9300 International and Maritime Telephone Number +1 (703) 527-3887

| Country | Organisation/Company | Address | Emergency number |
|----------------|---|---|-------------------------------------|
| United Kingdom | National Poisons Information Service (Newcastle Unit) | Claremont Place Newcastle-upon-Tyne, Newcastle | +44 191 2606182/+44 191 2606180 24H |

SECTION 2: HAZARDS IDENTIFICATION

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Repr. 1B H360Df

Aquatic Chronic 2 H411

Full text of hazard classes and H-statements : see section 16

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :



GHS08

GHS09

Signal word (CLP) : Danger

Hazardous ingredients : BBP, benzyl butyl phthalate

Hazard statements (CLP) : H360Df - May damage the unborn child. Suspected of damaging fertility
H411 - Toxic to aquatic life with long lasting effects

Precautionary statements (CLP) : P201 - Obtain special instructions before use
P202 - Do not handle until all safety precautions have been read and understood
P273 - Avoid release to the environment
P280 - Wear eye protection, protective gloves, protective clothing
P308+P313 - IF exposed or concerned: Get medical advice/attention

P391 - Collect spillage
 P405 - Store locked up
 P501 - Dispose of contents/container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation

Extra phrases : The product is seen as a medical device and therefore not subject to labelling (EU-regulation 1907/2006, article 2, paragraph 6c)

2.3. Other hazards

Other hazards not contributing to the classification : None under normal conditions.

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

3.1. Substance

Not applicable

3.2. Mixture

| Name | Product identifier | % | Classification according to Regulation (EC) No. 1272/2008 [CLP] |
|---|--|------------|--|
| BBP, benzyl butyl phthalate substance listed as REACH Candidate (Benzyl butyl phthalate (BBP)) substance listed in REACH Annex XIV (Benzyl butyl phthalate (BBP)) | (CAS No) 85-68-7 (EC no) 201-622-7 (EC index no) 607-430-00-3 (REACH-no) 01-2119489376-23 | =>5 - <25 | Repr. 1B, H360Df Aquatic Acute 1, H400 Aquatic Chronic 1, H410 |
| octamethylcyclotetrasiloxane | (CAS No) 556-67-2 (EC no) 209-136-7 (EC index no) 014-018-00-1 (REACH-no) 01-2119529238-36 | =>0,5 - <5 | Repr. 2, H361f Aquatic Chronic 4, H413 |
| sulfur | (CAS No) 7704-34-9 (EC no) 231-722-6 (EC index no) 016-094-00-1 (REACH-no) 01-2119487295-27 | =>1 - <5 | Skin Irrit. 2, H315 |

Full text of H-statements: see section 16

SECTION 4: FIRST AID MEASURES

4.1. Description of first aid measures

First-aid measures general : If you feel unwell, seek medical advice (show the label where possible).
 First-aid measures after inhalation : Remove victim to fresh air and keep at rest in a position comfortable for breathing. Get medical advice/attention if you feel unwell.
 First-aid measures after skin contact : Gently wash with plenty of soap and water. If skin irritation occurs: Get medical advice/attention.
 First-aid measures after eye contact : IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Get medical advice/attention if you feel unwell.
 First-aid measures after ingestion : If swallowed, rinse mouth with water (only if the person is conscious). Get medical advice/attention.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/injuries : May damage fertility or the unborn child.

4.3. Indication of any immediate medical attention and special treatment needed

No specific measures identified.

SECTION 5: FIREFIGHTING MEASURES

5.1. Extinguishing media

Suitable extinguishing media : Use extinguishing media appropriate for surrounding fire. Foam, carbon dioxide (CO2) and powder.
 Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

Fire hazard : Non flammable.
 Hazardous decomposition products in case of fire : Carbon dioxide. Carbon monoxide. Sulphur oxides. Halogenated compounds. metallic oxide.

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Avoid (reject) fire-fighting water to enter environment.

Protection during firefighting : Do not enter fire area without proper personal protective equipment, including respiratory protection.

SECTION 6: ACCIDENTAL RELEASE MEASURES

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Avoid contact with skin and eyes. Do not breathe vapour. Use personal protective equipment as required. No open flames. No smoking. Do not eat, drink or smoke in areas where product is used. Assure fresh air breathing.

For non-emergency personnel

Protective equipment : See Heading 8.

For emergency responders

No additional information available

6.2. Environmental precautions

Discharging into rivers and drains is forbidden.

6.3. Methods and material for containment and cleaning up

For containment : Collect all waste in suitable and labelled containers and dispose according to local legislation.

Methods for cleaning up : Soak up spills with inert solids, such as clay or diatomaceous earth as soon as possible.

6.4. Reference to other sections

For further information refer to section 13.

SECTION 7: HANDLING AND STORAGE

7.1. Precautions for safe handling

Precautions for safe handling : Do not eat, drink or smoke when using this product. Avoid contact with skin, eyes and clothing. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

Hygiene measures : Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Remove contaminated clothes. Wash contaminated clothing before reuse.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in a dry place. Store in a closed container. Store in original container. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

Storage area : Store in a well-ventilated place.

7.3. Specific end use(s)

Consult the supplier for further information.

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1. Control parameters

| BBP, benzyl butyl phthalate (85-68-7) | | |
|---------------------------------------|------------------------------|------------------------|
| United Kingdom | Local name | Benzyl butyl phthalate |
| United Kingdom | WEL TWA (mg/m ³) | 5 mg/m ³ |

8.2. Exposure controls

Appropriate engineering controls : Ensure good ventilation of the work station.

Personal protective equipment : Gloves. Safety glasses.

Hand protection : Wear suitable gloves. Nitrile rubber gloves. Layer thickness : 0,09mm. Breakthrough time : >480 min. STANDARD EN 374

Eye protection : Use splash goggles when eye contact due to splashing is possible. STANDARD EN 166

Skin and body protection : Wear suitable protective clothing

Respiratory protection : No special respiratory protection equipment is recommended under normal conditions of use with adequate ventilation



Other information : Do not eat, drink or smoke during use. Personal protective equipment should be chosen according to the CEN standards and in discussion with the supplier of the protective equipment. Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. For even minor contact, immediately remove contaminated clothing. Wash skin thoroughly with mild soap and water.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

9.1. Information on basic physical and chemical properties

Physical state : Liquid
 Appearance : Paste.
 Colour : white.
 Odour : Sulfur.
 Odour threshold : No data available
 pH : No data available
 Relative evaporation rate (butylacetate=1) : No data available
 Melting point : No data available
 Freezing point : No data available
 Boiling point : No data available
 Flash point : No data available
 Auto-ignition temperature : No data available
 Decomposition temperature : No data available
 Flammability (solid, gas) : No data available
 Vapour pressure : No data available
 Relative vapour density at 20 °C : No data available
 Relative density : > 1
 Solubility : Material insoluble in water.
 Log Pow : No data available
 Viscosity, kinematic : No data available
 Viscosity, dynamic : No data available
 Explosive properties : Product is not explosive.
 Oxidising properties : Non flammable.
 Explosive limits : No data available

9.2. Other information

Additional information : None to our knowledge.

SECTION 10: STABILITY AND REACTIVITY

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No polymerization.

10.4. Conditions to avoid

No flames. Eliminate all sources of ignition.

10.5. Incompatible materials

reducing materials.

10.6. Hazardous decomposition products

No decomposition if stored normally.

SECTION 11: TOXICOLOGICAL INFORMATION

11.1. Information on toxicological effects

Acute toxicity : Not classified

BBP, benzyl butyl phthalate (85-68-7)

| | |
|----------------------------|---------------|
| LD50 oral rat | 2330 mg/kg |
| LD50 dermal rat | 6700 mg/kg |
| LD50 dermal rabbit | > 10000 mg/kg |
| LC50 inhalation rat (mg/l) | > 6.7 mg/l/4h |

| octamethylcyclotetrasiloxane (556-67-2) | |
|--|---|
| LD50 oral rat | 1540 mg/kg |
| LD50 dermal rat | 1770 mg/kg |
| LC50 inhalation rat (mg/l) | 36 mg/l/4h |
| sulfur (7704-34-9) | |
| LD50 oral rat | > 3000 mg/kg |
| LD50 dermal rabbit | > 2000 mg/kg |
| LC50 inhalation rat (mg/l) | > 0.067 mg/l/4h |
| Skin corrosion/irritation | : Not classified |
| Serious eye damage/irritation | : Not classified |
| Respiratory or skin sensitisation | : Not classified |
| Germ cell mutagenicity | : Not classified |
| Carcinogenicity | : Not classified |
| Reproductive toxicity | : May damage the unborn child. Suspected of damaging fertility. |
| Specific target organ toxicity (single exposure) | : Not classified |
| Specific target organ toxicity (repeated exposure) | : Not classified |
| Aspiration hazard | : Not classified |

SECTION 12: ECOLOGICAL INFORMATION

12.1. Toxicity

Ecology - water : Toxic to aquatic life with long lasting effects.

| BBP, benzyl butyl phthalate (85-68-7) | |
|--|---|
| LC50 fish 1 | 0.51 mg/l (96 hours - Rainbow trout) |
| EC50 Daphnia 1 | 0.76 mg/l EC50 48h - Daphnia magna [mg/l] |
| EC50 other aquatic organisms 1 | 0.9 mg/l EC50 48h crustacea |
| IC50 algae | 0.1 mg/l (96 hours - Pseudokirchneriella subcapitata) |
| NOEC (chronic) | 0.06 mg/l |
| octamethylcyclotetrasiloxane (556-67-2) | |
| LC50 fish 1 | > 500 mg/l (96 hours - Brachydanio rerio, zebra-fish) |
| EC50 Daphnia 1 | 25.2 mg/l EC50 48h - Daphnia magna [mg/l] |
| sulfur (7704-34-9) | |
| LC50 fish 1 | 866 mg/l (96 hours - Brachydanio rerio, zebra-fish) |
| EC50 Daphnia 1 | > 5000 mg/l (48 hours - Daphnia magna) |
| IC50 algae | 12 mg/l |

12.2. Persistence and degradability

| Permlastic Base (Regular Body, Light Body & Heavy Body) | |
|--|-------------------------|
| Persistence and degradability | No data available. |
| BBP, benzyl butyl phthalate (85-68-7) | |
| Biodegradation | 93 % (OECD 301B method) |

12.3. Bioaccumulative potential

| Permlastic Base (Regular Body, Light Body & Heavy Body) | |
|--|----------|
| Bioaccumulative potential | No data. |
| BBP, benzyl butyl phthalate (85-68-7) | |
| Bioconcentration factor (BCF REACH) | 663 |
| Log Pow | 4.84 |
| octamethylcyclotetrasiloxane (556-67-2) | |
| Bioconcentration factor (BCF REACH) | 1698.24 |
| Log Pow | 5.09 |

12.4. Mobility in soil

| Permlastic Base (Regular Body, Light Body & Heavy Body) | |
|--|------------|
| Ecology - soil | Insoluble. |

12.5. Results of PBT and vPvB assessment

| Permlastic Base (Regular Body, Light Body & Heavy Body) | |
|--|--|
| This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII | |
| This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII | |

| Component | |
|---------------------------------------|---|
| BBP, benzyl butyl phthalate (85-68-7) | This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII |

12.6. Other adverse effects

Other adverse effects : None to our knowledge.
Additional information : No other effects known

SECTION 13: DISPOSAL CONSIDERATIONS

13.1. Waste treatment methods

Regional legislation (waste) : Dispose as hazardous waste.
Waste treatment methods : Recover the product with absorbent material. Dispose of contents/container in accordance with licensed collector's sorting instructions.
Waste disposal recommendations : Dispose in a safe manner in accordance with local/national regulations.
Ecology - waste materials : Avoid release to the environment.
European List of Waste (LoW) code : 18 01 06* - chemicals consisting of or containing dangerous substances

SECTION 14: TRANSPORT INFORMATION

In accordance with ADR / RID / IMDG / IATA / ADN

14.1. UN number

UN-No. (ADR) : 3077
UN-No. (IMDG) : 3077
UN-No. (IATA) : 3077
UN-No. (ADN) : 3077
UN-No. (RID) : 3077


14.2. UN proper shipping name

Proper Shipping Name (ADR) : ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S.
Proper Shipping Name (IMDG) : ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S.
Proper Shipping Name (IATA) : Environmentally hazardous substance, solid, n.o.s.
Proper Shipping Name (ADN) : ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S.
Proper Shipping Name (RID) : ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S.
Transport document description (ADR) : UN 3077 ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S. (zinc oxide(1314-13-2) ; BBP, benzyl butyl phthalate), 9, III, (E)
Transport document description (IMDG) : UN 3077 ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S., 9, III, MARINE POLLUTANT/ENVIRONMENTALLY HAZARDOUS

14.3. Transport hazard class(es)


ADR

Transport hazard class(es) (ADR) : 9
Danger labels (ADR) : 9




IMDG

Transport hazard class(es) (IMDG) : 9
Danger labels (IMDG) : 9



IATA


Transport hazard class(es) (IATA) : 9
Hazard labels (IATA) : 9



ADN

Transport hazard class(es) (ADN) : 9


Danger labels (ADN) : 9



RID

Transport hazard class(es) (RID) : 9

Danger labels (RID) : 9



14.4. Packing group

Packing group (ADR) : III

Packing group (IMDG) : III

Packing group (IATA) : III

Packing group (ADN) : III

Packing group (RID) : III

14.5. Environmental hazards

Dangerous for the environment : Yes

Marine pollutant : Yes

Other information : No supplementary information available

14.6. Special precautions for user

- Overland transport

- Classification code (ADR) : M7
- Special provisions (ADR) : 274, 335, 601, 375
- Limited quantities (ADR) : 5kg
- Excepted quantities (ADR) : E1
- Packing instructions (ADR) : P002, IBC08, LP02, R001
- Special packing provisions (ADR) : PP12, B3
- Mixed packing provisions (ADR) : MP10
- Portable tank and bulk container instructions (ADR) : T1, BK1, BK2
- Portable tank and bulk container special provisions (ADR) : TP33
- Tank code (ADR) : SGAV, LGBV
- Vehicle for tank carriage : AT
- Transport category (ADR) : 3
- Special provisions for carriage - Packages (ADR) : V13
- Special provisions for carriage - Bulk (ADR) : VC1, VC2
- Special provisions for carriage - Loading, unloading and handling (ADR) : CV13
- Hazard identification number (Kemler No.) : 90

Orange plates : 

Tunnel restriction code (ADR) : E
 EAC code : 2Z

- Transport by sea

Special provisions (IMDG) : 274, 335, 966, 967, 969
 Limited quantities (IMDG) : 5 kg
 Excepted quantities (IMDG) : E1
 Packing instructions (IMDG) : P002, LP02
 Special packing provisions (IMDG) : PP12
 IBC packing instructions (IMDG) : IBC08
 IBC special provisions (IMDG) : B3
 Tank instructions (IMDG) : T1, BK1, BK2, BK3
 Tank special provisions (IMDG) : TP33
 EmS-No. (Fire) : F-A
 EmS-No. (Spillage) : S-F
 Stowage category (IMDG) : A

- Air transport

PCA Excepted quantities (IATA) : E1
 PCA Limited quantities (IATA) : Y956
 PCA limited quantity max net quantity (IATA) : 30kgG
 PCA packing instructions (IATA) : 956
 PCA max net quantity (IATA) : 400kg
 CAO packing instructions (IATA) : 956
 CAO max net quantity (IATA) : 400kg
 Special provisions (IATA) : A97, A158, A179, A197
 ERG code (IATA) : 9L

Rail transport

No data available

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

SECTION 15: REGULATORY INFORMATION

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains a substance on the REACH candidate list in concentration $\geq 0.1\%$ or with a lower specific limit: Benzyl butyl phthalate (BBP) (EC 201-622-7, CAS 85-68-7)

Contains REACH Annex XIV substances:

| Substance name | Authorisation number | Sunset date | REACH authorisation exemptions |
|--|----------------------|-------------|--------------------------------|
| Benzyl butyl phthalate (BBP) (EC 201-622-7, CAS 85-68-7) | | 21/01/2015 | |

National regulations

EC-regulation 453/2010/EC, 1907/2006/EC (REACH), 1272/2008/EC (CLP), 790/2009/EC. Transport of dangerous goods (ADR/RID, IMDG, IATA/ICAO). Workplace exposure limits

15.2. Chemical safety assessment

No chemical safety assessment has been carried out for the substance or the mixture by the supplier

SECTION 16: OTHER INFORMATION

Indication of changes:

Hazards identification. Composition/information on ingredients.

| 2.1 | Classification according to Regulation (EC) No. | Modified | |
|-----|---|----------|--|
| | | | |

| | | | |
|-----|--|----------|--|
| | 1272/2008 [CLP] | | |
| 2.2 | Precautionary statements (CLP) | Modified | |
| 2.2 | Hazard statements (CLP) | Modified | |
| 3 | Composition/information on ingredients | Modified | |
| 5.2 | Hazardous decomposition products in case of fire | Modified | |

Date of issue : 24/02/2004
 Revision date : 04/03/2016
 Supersedes : 14/09/2015
 Version : 6.0
 Signature : A. Åsebø Murel

Full text of H- and EUH-statements:

| | |
|-------------------|---|
| Aquatic Acute 1 | Hazardous to the aquatic environment — Acute Hazard, Category 1 |
| Aquatic Chronic 1 | Hazardous to the aquatic environment — Chronic Hazard, Category 1 |
| Aquatic Chronic 2 | Hazardous to the aquatic environment — Chronic Hazard, Category 2 |
| Aquatic Chronic 4 | Hazardous to the aquatic environment — Chronic Hazard, Category 4 |
| Repr. 1B | Reproductive toxicity, Category 1B |
| Repr. 2 | Reproductive toxicity, Category 2 |
| Skin Irrit. 2 | Skin corrosion/irritation, Category 2 |
| H315 | Causes skin irritation |
| H360Df | May damage the unborn child. Suspected of damaging fertility |
| H361f | Suspected of damaging fertility |
| H400 | Very toxic to aquatic life |
| H410 | Very toxic to aquatic life with long lasting effects |
| H411 | Toxic to aquatic life with long lasting effects |
| H413 | May cause long lasting harmful effects to aquatic life |

The information in this safety data sheet is based on information from the manufacturer/supplier, present European and national legislation, and presupposes that the product is used within the specified area of application.



Safety Data Sheet Permlastic Catalyst (Light Body)

SECTION 1: IDENTIFICATION OF THE SUBSTANCE/MIXTURE AND OF THE COMPANY/UNDERTAKING

1.1. Product identifier

Product name : Permlastic Catalyst (Light Body)

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses

Main use category : Professional use

Use of the substance/mixture : Preparation intended for dental medical use

Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

KERRHAWE S.A.
Via Strecce n°4
6934 Bioggio (Switzerland)
T 00-800-41-050-505

Manufacturer

Kerr Italia S.r.l.
Via Passanti, 332
84018 Scafati (SA) - Italy
T +39-081-850-8311

Contact person : safety@kerrhawe.com - tel. 00-800-41-050-505 (08.00-17.00)

1.4. Emergency telephone number

Emergency number : CHEMTREC® Emergency Call Center. Emergency Telephone Number (for USA only) 001-800-424-9300 International and Maritime Telephone Number +1 (703) 527-3887

| Country | Organisation/Company | Address | Emergency number |
|----------------|--|---|--|
| Gibraltar | GHA Call Centre Zone 2, Level3, St Bernard's Hospital | Harbour Views Road | +350 200 79700 +350 200 72266 |
| Ireland | National Poisons Information Centre Beaumont Hospital | PO Box 1297 Beaumont Road 9 Dublin | +353 1 809 2566 (Healthcare professionals- 24/7) +353 1 809 2166 (public, 8am - 10pm, 7/7) |
| Malta | Medicines & Poisons Info Office | Mater Dei Hospital MSD Msida | +356 2545 6504 |
| United Kingdom | National Poisons Information Service (Newcastle Unit) | Claremont Place Newcastle-upon-Tyne, Newcastle | +44 191 2606182/+44 191 2606180 24H |

SECTION 2: HAZARDS IDENTIFICATION

2.1. Classification of the substance or mixture



Classification according to Regulation (EC) No. 1272/2008 [CLP]

Repr. 1A H360
STOT RE 2 H373
Aquatic Acute 1 H400
Aquatic Chronic 1 H410

Full text of H statements : see section 16

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

| | | |
|--------------------------------|---|---|
| Hazard pictograms (CLP) | : |   |
| | | GHS08 GHS09 |
| Signal word (CLP) | : | Danger |
| Hazardous ingredients | : | lead dioxide |
| Hazard statements (CLP) | : | H360 - May damage fertility or the unborn child. H373 - May cause damage to organs through prolonged or repeated exposure. H410 - Very toxic to aquatic life with long lasting effects. |
| Precautionary statements (CLP) | : | P201 - Obtain special instructions before use. P202 - Do not handle until all safety precautions have been read and understood. P260 - Do not breathe vapours, dust. P273 - Avoid release to the environment. P280 - Wear eye protection, protective gloves, protective clothing. P314 - Get medical advice/attention if you feel unwell. P391 - Collect spillage. P405 - Store locked up. P501 - Dispose of contents and container to hazardous or special waste collection point, in accordance with local, regional, national and/or international regulation. |
| Extra phrases | : | The product is seen as a medical device and therefore not subject to labelling (EU-regulation 1272/2008, article 1, paragraph 5d). A safety data sheet is not required for this product under Article 31 of REACH. This Product Safety Information Sheet has been created on a voluntary basis |

2.3. Other hazards

Other hazards not contributing to the classification : None under normal conditions.

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

3.1. Substances

Not applicable

3.2. Mixtures

| Name | Product identifier | % | Classification according to Regulation (EC) No. 1272/2008 [CLP] |
|--------------|---|---------|--|
| zinc oxide | (CAS-No.) 1314-13-2 (EC-No.) 215-222-5 (EC Index-No.) 030-013-00-7 (REACH-no) 01-2119463881-32 | 10 - 50 | Aquatic Acute 1, H400 (M=10) Aquatic Chronic 1, H410 |
| lead dioxide | (CAS-No.) 1309-60-0 (EC-No.) 215-174-5 (REACH-no) 01-2119958814-25 | 5 - 25 | Ox. Sol. 2, H272 Acute Tox. 4 (Oral), H302 Acute Tox. 4 (Inhalation:dust,mist), H332 Repr. 1A, H360 STOT RE 2, H373 Aquatic Acute 1, H400 Aquatic Chronic 1, H410 |

Full text of H-statements: see section 16

SECTION 4: FIRST AID MEASURES

4.1. Description of first aid measures

| | | |
|---------------------------------------|---|---|
| First-aid measures general | : | IF exposed or concerned: Get medical advice/attention. |
| First-aid measures after inhalation | : | Remove person to fresh air and keep comfortable for breathing. Get medical advice/attention if you feel unwell. |
| First-aid measures after skin contact | : | Gently wash with plenty of soap and water. If skin irritation occurs: Get medical advice/attention. |
| First-aid measures after eye contact | : | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Get medical advice/attention if you feel unwell. |

First-aid measures after ingestion : If swallowed, rinse mouth with water (only if the person is conscious). Get medical advice/attention.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects : May cause damage to organs. May damage fertility. May damage the unborn child.

4.3. Indication of any immediate medical attention and special treatment needed

No specific measures identified. In all cases of doubt, or when symptoms persist, seek medical attention.

SECTION 5: FIREFIGHTING MEASURES

5.1. Extinguishing media

Suitable extinguishing media : Use extinguishing media appropriate for surrounding fire. Foam, carbon dioxide (CO₂) and powder.

Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

Fire hazard : Non flammable.

Hazardous decomposition products in case of fire : Carbon dioxide. Carbon monoxide. metallic oxide. Nitrogen oxides. Lead compounds :

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Prevent fire fighting water from entering the environment.

Protection during firefighting : Do not enter fire area without proper personal protective equipment, including respiratory protection.

SECTION 6: ACCIDENTAL RELEASE MEASURES

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Avoid contact with skin and eyes. Use personal protective equipment as required. Avoid all unnecessary exposure. Do not eat, drink or smoke in areas where product is used. Ensure adequate ventilation, especially in confined areas.

For non-emergency personnel

Protective equipment : See Heading 8.

Emergency procedures : Evacuate unnecessary personnel.

For emergency responders

No additional information available

6.2. Environmental precautions

Discharging into rivers and drains is forbidden. Notify authorities if liquid enters sewers or public waters.

6.3. Methods and material for containment and cleaning up

For containment : Collect all waste in suitable and labelled containers and dispose according to local legislation.

Methods for cleaning up : Soak up spills with inert solids, such as clay or diatomaceous earth as soon as possible.

6.4. Reference to other sections

For further information refer to section 13. See Heading 8. Exposure controls and personal protection.

SECTION 7: HANDLING AND STORAGE

7.1. Precautions for safe handling

Precautions for safe handling : Do not eat, drink or smoke when using this product. Avoid contact with skin, eyes and clothing. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood.

Hygiene measures : Avoid all eye and skin contact and do not breathe vapour and mist. Always wash hands after handling the product. Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Wash contaminated clothing before reuse. Contaminated work clothing should not be allowed out of the workplace.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in a dry place. Store in a closed container. Store in original container. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.

Incompatible materials : reducing materials. Food.

Storage area : Store in a well-ventilated place.

7.3. Specific end use(s)

For professional use only.

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1. Control parameters

| lead dioxide (1309-60-0) | | |
|--------------------------|---------------------------|----------------------------------|
| EU | Local name | Lead and its inorganic compounds |
| EU | IOELV TWA (mg/m³) | 100 µg/m³ |
| EU | Notes | SCOEL Recommendations (2002) |
| Ireland | Local name | Lead |
| Ireland | OEL (8 hours ref) (mg/m³) | 0.15 mg/m³ |
| Ireland | Notes (IE) | Repr.1A, BOELV |
| zinc oxide (1314-13-2) | | |
| EU | Local name | Zinc oxide |
| EU | Notes | (Ongoing) |
| Ireland | Local name | Zinc oxide, fume |
| Ireland | OEL (8 hours ref) (mg/m³) | 2 mg/m³ R (Respirable) |
| Ireland | OEL (15 min ref) (mg/m3) | 10 mg/m³ |

8.2. Exposure controls

- Appropriate engineering controls : Ensure good ventilation of the work station.
- Personal protective equipment : Gloves.
- Hand protection : Wear suitable gloves. Nitrile rubber gloves. Layer thickness : 0,09mm. Breakthrough time : >480 min. STANDARD EN 374.
- Eye protection : Use splash goggles when eye contact due to splashing is possible. STANDARD EN 166.
- Skin and body protection : Wear suitable protective clothing
- Respiratory protection : No special respiratory protection equipment is recommended under normal conditions of use with adequate ventilation. If dust are formed : Extra personal protection: P3 filter respirator for toxic particles. Standard EN 143.



- Other information : Do not eat, drink or smoke during use. Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

9.1. Information on basic physical and chemical properties

- Physical state : Solid
- Appearance : Paste.
- Colour : brown. Violet. Grey.
- Odour : Fruity.
- Odour threshold : Not determined
- pH : Not determined
- Relative evaporation rate (butylacetate=1) : Not determined
- Melting point : Not determined
- Freezing point : Not determined
- Boiling point : Not determined
- Flash point : Not determined
- Auto-ignition temperature : Not determined
- Decomposition temperature : Not determined
- Flammability (solid, gas) : Non flammable.
- Vapour pressure : Not determined
- Relative vapour density at 20 °C : Not determined
- Relative density : > 1
- Solubility : Material insoluble in water.
- Log Pow : Not determined
- Viscosity, kinematic : Not determined

| | |
|----------------------|-----------------------------|
| Viscosity, dynamic | : Not determined |
| Explosive properties | : Product is not explosive. |
| Oxidising properties | : Non flammable. |
| Explosive limits | : Not determined |

9.2. Other information

Additional information : None to our knowledge.

SECTION 10: STABILITY AND REACTIVITY**10.1. Reactivity**

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No polymerization.

10.4. Conditions to avoid

No flames. Eliminate all sources of ignition. Keep away from heat and direct sunlight.

10.5. Incompatible materials

reducing materials. oxidizing materials.

10.6. Hazardous decomposition products

No decomposition if stored and used normally.

SECTION 11: TOXICOLOGICAL INFORMATION**11.1. Information on toxicological effects**

Acute toxicity : Not classified

| | |
|-------------------------------|-------------------------|
| ATE CLP (oral) | 2512.6 mg/kg bodyweight |
| zinc oxide (1314-13-2) | |
| LD50 oral rat | > 15000 mg/kg |
| LD50 dermal rat | > 2000 mg/kg |
| LC50 inhalation rat (mg/l) | 5.7 mg/l/4h |

| | |
|-----------------------------------|--|
| Skin corrosion/irritation | : Not classified Based on available data, the classification criteria are not met pH: Not determined |
| Serious eye damage/irritation | : Not classified Based on available data, the classification criteria are not met pH: Not determined |
| Respiratory or skin sensitisation | : Not classified Based on available data, the classification criteria are not met |
| Germ cell mutagenicity | : Not classified Based on available data, the classification criteria are not met |
| Carcinogenicity | : Not classified Based on available data, the classification criteria are not met |
| Reproductive toxicity | : May damage fertility or the unborn child. |
| STOT-single exposure | : Not classified Based on available data, the classification criteria are not met |
| STOT-repeated exposure | : May cause damage to organs through prolonged or repeated exposure. |
| Aspiration hazard | : Not classified Based on available data, the classification criteria are not met |

SECTION 12: ECOLOGICAL INFORMATION**12.1. Toxicity**

Ecology - general : Very toxic to aquatic life with long lasting effects.

| | |
|-------------------------------|--|
| zinc oxide (1314-13-2) | |
| LC50 fish 1 | 1.1 mg/l (96 hours - Rainbow trout) |
| EC50 Daphnia 1 | 98 µg/l (48 hours - Daphnia magna) |
| NOEC (chronic) | 0.017 mg/l 72 hours - Pseudokirchnerella subcapitata |

12.2. Persistence and degradability

| Permlastic Catalyst (Light Body) | |
|----------------------------------|--------------------|
| Persistence and degradability | No data available. |

12.3. Bioaccumulative potential

| Permlastic Catalyst (Light Body) | |
|-------------------------------------|----------------|
| Log Pow | Not determined |
| Bioaccumulative potential | No data. |
| lead dioxide (1309-60-0) | |
| Bioconcentration factor (BCF REACH) | 70000 |
| zinc oxide (1314-13-2) | |
| Bioconcentration factor (BCF REACH) | 60960 |

12.4. Mobility in soil

| Permlastic Catalyst (Light Body) | |
|----------------------------------|------------|
| Ecology - soil | Insoluble. |

12.5. Results of PBT and vPvB assessment

| Permlastic Catalyst (Light Body) | |
|--|--|
| This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII | |
| This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII | |

12.6. Other adverse effects

Other adverse effects : None to our knowledge.
 Additional information : Avoid release to the environment.

SECTION 13: DISPOSAL CONSIDERATIONS

13.1. Waste treatment methods

Regional legislation (waste) : Dispose as hazardous waste.
 Waste treatment methods : Dispose of contents/container in accordance with licensed collector's sorting instructions. Do not discharge into drains.
 Product/Packaging disposal recommendations : Dispose in a safe manner in accordance with local/national regulations. Dispose of contents/container to a hazardous or special waste collection point.
 Ecology - waste materials : Avoid release to the environment.
 European List of Waste (LoW) code : 18 01 06* - chemicals consisting of or containing dangerous substances

SECTION 14: TRANSPORT INFORMATION

In accordance with ADR / RID / IMDG / IATA / ADN


| ADR | IMDG | IATA | ADN | RID |
|--|---|--|--|--|
| 14.1. UN number | | | | |
| 3077 | 3077 | 3077 | 3077 | 3077 |
| 14.2. UN proper shipping name | | | | |
| ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S. | ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S. | Environmentally hazardous substance, solid, n.o.s. | ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S. | ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S. |
| Transport document description | | | | |
| UN 3077 ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S. (zinc oxide), 9, III, (-) | UN 3077 ENVIRONMENTALLY HAZARDOUS SUBSTANCE, SOLID, N.O.S. (zinc oxide), 9, III, MARINE POLLUTANT | | | |
| 14.3. Transport hazard class(es) | | | | |
| 9 | 9 | 9 | 9 | 9 |
| | | | | |

| ADR | IMDG | IATA | ADN | RID |
|--|---|-------------------------------------|-------------------------------------|-------------------------------------|
| 14.4. Packing group | | | | |
| III | III | III | III | III |
| 14.5. Environmental hazards | | | | |
| Dangerous for the environment : Yes | Dangerous for the environment : Yes Marine pollutant : Yes | Dangerous for the environment : Yes | Dangerous for the environment : Yes | Dangerous for the environment : Yes |
| No supplementary information available | | | | |

14.6. Special precautions for user

- Overland transport

Special provisions (ADR) : 274, 335, 375, 601
 Limited quantities (ADR) : 5kg
 Excepted quantities (ADR) : E1
 Hazard identification number (Kemler No.) : 90
 Orange plates :



EAC code : 2Z

- Transport by sea

Special provisions (IMDG) : 274, 335, 966, 967, 969
 EmS-No. (Fire) : F-A
 EmS-No. (Spillage) : S-F

- Air transport

PCA Excepted quantities (IATA) : E1
 PCA Limited quantities (IATA) : Y956
 Special provisions (IATA) : A97, A158, A179, A197

Rail transport

No data available

14.7. Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

SECTION 15: REGULATORY INFORMATION

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

EU-Regulations

Contains no REACH substances with Annex XVII restrictions
 Contains no substance on the REACH candidate list
 Contains no REACH Annex XIV substances

National regulations

EC-regulation 2015/830 /EC, 1907/2006/EC (REACH), 1272/2008/EC (CLP), 790/2009/EC. Transport of dangerous goods (ADR/RID, IMDG, IATA/ICAO). Workplace exposure limits.

15.2. Chemical safety assessment

No chemical safety assessment has been carried out for the substance or the mixture by the supplier

SECTION 16: OTHER INFORMATION

Indication of changes:

| | | | |
|------|---|----------|--|
| 2.1 | Classification according to Regulation (EC) No. 1272/2008 [CLP] | Modified | |
| 2.2 | Precautionary statements (CLP) | Modified | |
| 2.2 | Hazard statements (CLP) | Modified | |
| 3 | Composition/information on ingredients | Modified | |
| 10.5 | Incompatible materials | Modified | |

Date of issue : 03/03/2004
 Revision date : 30/12/2017
 Supersedes : 03/01/2017
 Date of total revision : 30/12/2017
 Version : 6.0
 Signature : A. Åsebø Murel

Full text of H- and EUH-statements:

| | |
|-------------------------------------|--|
| Acute Tox. 4 (Inhalation:dust,mist) | Acute toxicity (inhalation:dust,mist) Category 4 |
| Acute Tox. 4 (Oral) | Acute toxicity (oral), Category 4 |
| Aquatic Acute 1 | Hazardous to the aquatic environment — Acute Hazard, Category 1 |
| Aquatic Chronic 1 | Hazardous to the aquatic environment — Chronic Hazard, Category 1 |
| Ox. Sol. 2 | Oxidising Solids, Category 2 |
| Repr. 1A | Reproductive toxicity, Category 1A |
| STOT RE 2 | Specific target organ toxicity — Repeated exposure, Category 2 |
| H272 | May intensify fire; oxidiser. |
| H302 | Harmful if swallowed. |
| H332 | Harmful if inhaled. |
| H360 | May damage fertility or the unborn child. |
| H373 | May cause damage to organs through prolonged or repeated exposure. |
| H400 | Very toxic to aquatic life. |
| H410 | Very toxic to aquatic life with long lasting effects. |

The information in this safety data sheet is based on information from the manufacturer/supplier, present european and national legislation, and presupposes that the product is used within the specified area of application.