

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Occlusion spray non flammable green, red, blue**

Revision date: 30.05.2022

Page 1 of 13

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1. Product identifier**

Occlusion spray non flammable green, red, blue

**Product code:**

11140-1

11145-1

11146-1

UFI: WSU2-H0FH-P00D-T6GE

**1.2. Relevant identified uses of the substance or mixture and uses advised against****1.3. Details of the supplier of the safety data sheet**

Company name:	DFS-DIAMON GmbH	
Street:	Ländenstr. 1	
Place:	D-93339 Riedenburg	
Telephone:	+49 (0)9442 9189 0	Telefax: +49 (0)9442 9189 37
e-mail:	info@dfs-diamon.de	
Contact person:	Armin Aunkofer	
e-mail:	armin.aunkofer@dfs-diamon.de	
Internet:	www.dfs-diamon.de	

**1.4. Emergency telephone number:** 24 hour Emergency Contact: +49 (0) 700 24112112 (DFS)  
+1 872 5888271

**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture****Regulation (EC) No 1272/2008**

Aerosol 3; H229  
STOT SE 3; H336  
Aquatic Chronic 3; H412

Full text of hazard statements: see SECTION 16.

**2.2. Label elements****Regulation (EC) No 1272/2008****Hazard components for labelling**

pentane

**Signal word:** Warning**Pictograms:****Hazard statements**

H229	Pressurised container: May burst if heated.
H336	May cause drowsiness or dizziness.
H412	Harmful to aquatic life with long lasting effects.

**Precautionary statements**

P210	Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P251	Do not pierce or burn, even after use.
P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P312	Call a POISON CENTER/doctor if you feel unwell.

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## Occlusion spray non flammable green, red, blue

Revision date: 30.05.2022

Page 2 of 13

P403+P233 Store in a well-ventilated place. Keep container tightly closed.  
 P410+P412 Protect from sunlight. Do not expose to temperatures exceeding 50 °C/122 °F.  
 P501 Dispose of waste according to applicable legislation

### Special labelling of certain mixtures

10 % by mass of the contents are flammable.

### 2.3. Other hazards

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures

#### Hazardous components

CAS No	Chemical name			Quantity
	EC No	Index No	REACH No	
	Classification (Regulation (EC) No 1272/2008)			
64-17-5	ethanol, ethyl alcohol			1 - < 5 %
	200-578-6	603-002-00-5	01-2119457610-43	
	Flam. Liq. 2, Eye Irrit. 2; H225 H319			
109-66-0	pentane			1 - < 5 %
	203-692-4	601-006-00-1	01-2119459286-30	
	Flam. Liq. 2, STOT SE 3, Asp. Tox. 1, Aquatic Chronic 2; H225 H336 H304 H411 EUH066			
78-93-3	butanone; ethyl methyl ketone			< 0.1 %
	201-159-0	606-002-00-3	01-2119457290-43	
	Flam. Liq. 2, Eye Irrit. 2, STOT SE 3; H225 H319 H336 EUH066			

Full text of H and EUH statements: see section 16.

#### Specific Conc. Limits, M-factors and ATE

CAS No	EC No	Chemical name	Quantity
	Specific Conc. Limits, M-factors and ATE		
64-17-5	200-578-6	ethanol, ethyl alcohol	1 - < 5 %
	inhalation: LC50 = 117 mg/l (vapours); oral: LD50 = 10470 mg/kg Eye Irrit. 2; H319: >= 50 - 100		
78-93-3	201-159-0	butanone; ethyl methyl ketone	< 0.1 %
	dermal: LD50 = 6480 mg/kg; oral: LD50 = 2193 mg/kg		

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### General information

When in doubt or if symptoms are observed, get medical advice.

#### After inhalation

Provide fresh air. In case of accident or unwellness, seek medical advice immediately (show directions for use or safety data sheet if possible).

#### After contact with skin

Wash with plenty of water. Take off contaminated clothing and wash it before reuse. In case of skin reactions, consult a physician.

#### After contact with eyes

Rinse immediately carefully and thoroughly with eye-bath or water. Remove contact lenses, if present and easy to do. Continue rinsing. In case of eye irritation consult an ophthalmologist.

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## Occlusion spray non flammable green, red, blue

Revision date: 30.05.2022

Page 3 of 13

### After ingestion

Rinse mouth immediately and drink plenty of water. Provide fresh air.

### 4.2. Most important symptoms and effects, both acute and delayed

No information available.

### 4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

#### Suitable extinguishing media

Co-ordinate fire-fighting measures to the fire surroundings.

### 5.2. Special hazards arising from the substance or mixture

Non-flammable. Vapours can form explosive mixtures with air. Pressurised container: May burst if heated.

In case of fire may be liberated: Pyrolysis products, containing fluorine, Pyrolysis products, toxic.

### 5.3. Advice for firefighters

In case of fire: Wear self-contained breathing apparatus. Full protective suit.

#### Additional information

Use water spray jet to protect personnel and to cool endangered containers. Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### General advice

Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Evacuate area.

#### For non-emergency personnel

Provide adequate ventilation. Use personal protection equipment.

#### For emergency responders

Use personal protection equipment.

### 6.2. Environmental precautions

Do not allow to enter into surface water or drains.

### 6.3. Methods and material for containment and cleaning up

#### For containment

Stop leak if safe to do so.

#### For cleaning up

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal. Ventilate affected area.

#### Other information

Clean contaminated articles and floor according to the environmental legislation.

### 6.4. Reference to other sections

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

#### Advice on safe handling

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Occlusion spray non flammable green, red, blue

Revision date: 30.05.2022

Page 4 of 13

clothes. Use personal protection equipment.

#### Advice on protection against fire and explosion

Heating causes rise in pressure with risk of bursting. Vapours can form explosive mixtures with air.  
Protect from sunlight. Do not expose to temperatures exceeding 50 °C/122 °F.

#### Advice on general occupational hygiene

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme.  
Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

#### Further information on handling

Do not pierce or burn, even after use.

#### 7.2. Conditions for safe storage, including any incompatibilities

##### Requirements for storage rooms and vessels

Keep container tightly closed. Keep in a cool, well-ventilated place.

##### Hints on joint storage

Do not store together with: Pyrophoric or self-heating substances.

##### Further information on storage conditions

Keep away from heat. Protect from sunlight.

### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

##### Occupational exposure limits

CAS No	Substance	ppm	mg/m <sup>3</sup>	fib/cm <sup>3</sup>	Category	Origin
-	Dust non-specific, respirable	-	4		TWA (8 h)	
-	Dust non-specific, total inhalable	-	10		TWA (8 h)	
64-17-5	Ethanol	1000	-		STEL (15 min)	
78-93-3	Methyl ethyl ketone (MEK) (Butan-2-one)	200	600		TWA (8 h)	
		300	900		STEL (15 min)	
109-66-0	n-Pentane	1000	-		TWA (8 h)	

##### Biological limit values

CAS No	Substance	Parameter	Value	Test material	Sampling time
78-93-3	Butan-2-one	Butan-2-one	70 µmol/L	Urine	Post shift

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Occlusion spray non flammable green, red, blue**

Revision date: 30.05.2022

Page 5 of 13

**DNEL/DMEL values**

CAS No	Substance			
DNEL type	Exposure route	Effect	Value	
64-17-5	ethanol, ethyl alcohol			
Worker DNEL, long-term	inhalation	systemic	950 mg/m <sup>3</sup>	
Worker DNEL, acute	inhalation	local	1900 mg/m <sup>3</sup>	
Worker DNEL, long-term	dermal	systemic	343 mg/kg bw/day	
Consumer DNEL, long-term	inhalation	systemic	114 mg/m <sup>3</sup>	
Consumer DNEL, acute	inhalation	local	950 mg/m <sup>3</sup>	
Consumer DNEL, long-term	dermal	systemic	206 mg/kg bw/day	
Consumer DNEL, acute	oral	systemic	87 mg/kg bw/day	
109-66-0	pentane			
Worker DNEL, long-term	inhalation	systemic	3000 mg/m <sup>3</sup>	
Worker DNEL, long-term	dermal	systemic	432 mg/kg bw/day	
Consumer DNEL, long-term	inhalation	systemic	643 mg/m <sup>3</sup>	
Consumer DNEL, long-term	dermal	systemic	214 mg/kg bw/day	
Consumer DNEL, long-term	oral	systemic	214 mg/kg bw/day	
78-93-3	butanone; ethyl methyl ketone			
Worker DNEL, long-term	inhalation	systemic	600 mg/m <sup>3</sup>	
Worker DNEL, long-term	dermal	systemic	1161 mg/kg bw/day	
Consumer DNEL, long-term	inhalation	systemic	106 mg/m <sup>3</sup>	
Consumer DNEL, long-term	dermal	systemic	412 mg/kg bw/day	
Consumer DNEL, long-term	oral	systemic	31 mg/kg bw/day	

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Occlusion spray non flammable green, red, blue

Revision date: 30.05.2022

Page 6 of 13

#### PNEC values

CAS No	Substance	Value
Environmental compartment		
64-17-5	ethanol, ethyl alcohol	
Freshwater		0,96 mg/l
Marine water		0,79 mg/l
Freshwater sediment		3,6 mg/kg
Marine sediment		2,9 mg/kg
Secondary poisoning		0,72 mg/kg
Micro-organisms in sewage treatment plants (STP)		580 mg/l
Soil		0,63 mg/kg
109-66-0	pentane	
Freshwater		0,23 mg/l
Freshwater (intermittent releases)		0,88 mg/l
Marine water		0,23 mg/l
Freshwater sediment		1,2 mg/kg
Marine sediment		1,2 mg/kg
Soil		0,55 mg/kg
78-93-3	butanone; ethyl methyl ketone	
Freshwater		55,8 mg/l
Marine water		55,8 mg/l
Micro-organisms in sewage treatment plants (STP)		709 mg/l
Freshwater sediment		284,74 mg/kg
Marine sediment		284,7 mg/kg
Soil		22,5 mg/kg
Secondary poisoning		1000 mg/kg

#### 8.2. Exposure controls



##### Appropriate engineering controls

Provide adequate ventilation as well as local exhaustion at critical locations.

##### Individual protection measures, such as personal protective equipment

###### Eye/face protection

Use eye protection according to EN 166.

###### Hand protection

Wear suitable gloves tested to EN374.

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

###### Skin protection

Wear suitable protective clothing.

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Occlusion spray non flammable green, red, blue**

Revision date: 30.05.2022

Page 7 of 13

**Respiratory protection**

In case of inadequate ventilation wear respiratory protection.

Respiratory protection necessary at: exceeding exposure limit values, aerosol or mist formation.

**Thermal hazards**

Heating causes rise in pressure with risk of bursting.

**Environmental exposure controls**

Avoid release to the environment.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state:	Liquid (Aerosol)
Colour:	According to product specification
Odour:	characteristic
Odour threshold:	not determined

**Changes in the physical state**

Melting point/freezing point:	not determined
Boiling point or initial boiling point and boiling range:	(Propellant) -19 °C
Flash point:	not determined (*)

**Flammability**

Solid/liquid:	not applicable
Gas:	not applicable
Lower explosion limits:	not determined
Upper explosion limits:	not determined
Auto-ignition temperature:	(Propellant) 368 °C
Decomposition temperature:	not determined
pH-Value:	not determined
Viscosity / kinematic:	not determined
Water solubility:	practically insoluble

**Solubility in other solvents**

not determined

Partition coefficient n-octanol/water:	not determined
Vapour pressure: (at 20 °C)	5740 hPa
Vapour pressure: (at 50 °C)	13200 hPa
Density:	not determined
Relative vapour density:	not determined
Particle characteristics:	not determined

**9.2. Other information****Other safety characteristics**

(\*) Propellant: flame-resistant

Danger of explosion: Heating causes rise in pressure with risk of bursting.

**SECTION 10: Stability and reactivity**

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Occlusion spray non flammable green, red, blue

Revision date: 30.05.2022

Page 8 of 13

**10.1. Reactivity**

Aerosol

**10.2. Chemical stability**

The product is stable under storage at normal ambient temperatures.

**10.3. Possibility of hazardous reactions**

Vapours can form explosive mixtures with air.  
Heating causes rise in pressure with risk of bursting.

**10.4. Conditions to avoid**

Protect from direct sunlight. Keep away from heat.

**10.5. Incompatible materials**

Pyrophoric or self-heating substances.

**10.6. Hazardous decomposition products**

In case of fire may be liberated: Pyrolysis products, containing fluorine, Pyrolysis products, toxic.

## SECTION 11: Toxicological information

**11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008****Acute toxicity**

Based on available data, the classification criteria are not met.

CAS No	Chemical name				
	Exposure route	Dose	Species	Source	Method
64-17-5	ethanol, ethyl alcohol				
	oral	LD50 mg/kg	10470	Rat	Manufacturer
	inhalation (4 h) vapour	LC50	117 mg/l	Rat	Manufacturer
78-93-3	butanone; ethyl methyl ketone				
	oral	LD50 mg/kg	2193	Rat	OECD 423
	dermal	LD50 mg/kg	6480	Rabbit	GESTIS

**Irritation and corrosivity**

Based on available data, the classification criteria are not met.

**Sensitising effects**

Based on available data, the classification criteria are not met.

**Carcinogenic/mutagenic/toxic effects for reproduction**

Based on available data, the classification criteria are not met.

**STOT-single exposure**

May cause drowsiness or dizziness. (pentane)

**STOT-repeated exposure**

Based on available data, the classification criteria are not met.

**Aspiration hazard**

Based on available data, the classification criteria are not met.

**Information on likely routes of exposure**

Eye contact, Skin contact, Inhalation.

Active agent: oral

**11.2. Information on other hazards****Endocrine disrupting properties**

No information available.



## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Occlusion spray non flammable green, red, blue

Revision date: 30.05.2022

Page 9 of 13

#### SECTION 12: Ecological information

##### 12.1. Toxicity

Harmful to aquatic life with long lasting effects.

CAS No	Chemical name					
	Aquatic toxicity	Dose	[h]   [d]	Species	Source	Method
64-17-5	ethanol, ethyl alcohol					
	Acute fish toxicity	LC50 mg/l	14200	96 h	Pimephales promelas (fathead minnow)	Manufacturer
	Acute crustacea toxicity	EC50 14221 mg/l	9268 -	48 h	Daphnia magna	Manufacturer
78-93-3	butanone; ethyl methyl ketone					
	Acute fish toxicity	LC50 mg/l	2993	96 h	Pimephales promelas (fathead minnow)	OECD 203
	Acute algae toxicity	ErC50 mg/l	2029	96 h	Selenastrum capricornutum	OECD 201
	Acute crustacea toxicity	EC50	308 mg/l	48 h	Daphnia magna (Big water flea)	OECD 202
	Fish toxicity	NOEC mg/l	1170	4 d	Pimephales promelas (fathead minnow)	OECD 203
	Algae toxicity	NOEC mg/l	1240	4 d	Selenastrum capricornutum	OECD 201
	Crustacea toxicity	NOEC	68 mg/l	2 d	Daphnia magna (Big water flea)	OECD 202

##### 12.2. Persistence and degradability

The product has not been tested.

CAS No	Chemical name			
	Method	Value	d	Source
	Evaluation			
64-17-5	ethanol, ethyl alcohol			
	OECD 301B/ ISO 9439/ EEC 92/69/V, C.4-C	97 %	28	IUCLID
	Readily biodegradable (according to OECD criteria).			
78-93-3	butanone; ethyl methyl ketone			
		98 %	28	OECD 301D/ EEC 92/69/V, C.4-E
	Readily biodegradable (according to OECD criteria).			

##### 12.3. Bioaccumulative potential

The product has not been tested.

##### Partition coefficient n-octanol/water

CAS No	Chemical name	Log Pow
64-17-5	ethanol, ethyl alcohol	-0,31
78-93-3	butanone; ethyl methyl ketone	0,3

##### 12.4. Mobility in soil

The product has not been tested.

##### 12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.

##### 12.6. Endocrine disrupting properties

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Occlusion spray non flammable green, red, blue

Revision date: 30.05.2022

Page 10 of 13

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

#### 12.7. Other adverse effects

No information available.

#### Further information

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

##### Disposal recommendations

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

##### Contaminated packaging

Do not pierce or burn, even after use. Handle contaminated packages in the same way as the substance itself.

### SECTION 14: Transport information

#### Land transport (ADR/RID)

<b>14.1. UN number or ID number:</b>	UN 1950
<b>14.2. UN proper shipping name:</b>	AEROSOLS
<b>14.3. Transport hazard class(es):</b>	2
<b>14.4. Packing group:</b>	-
Hazard label:	2.2



Classification code:	5A
Special Provisions:	190 327 344 625
Limited quantity:	1 L
Excepted quantity:	E0
Transport category:	3
Tunnel restriction code:	E

#### Inland waterways transport (ADN)

<b>14.1. UN number or ID number:</b>	UN 1950
<b>14.2. UN proper shipping name:</b>	AEROSOLS
<b>14.3. Transport hazard class(es):</b>	2
<b>14.4. Packing group:</b>	-
Hazard label:	2.2



Classification code:	5A
Special Provisions:	190 327 344 625
Limited quantity:	1 L
Excepted quantity:	E0

#### Marine transport (IMDG)

<b>14.1. UN number or ID number:</b>	UN 1950
<b>14.2. UN proper shipping name:</b>	AEROSOLS
<b>14.3. Transport hazard class(es):</b>	2.2
<b>14.4. Packing group:</b>	-
Hazard label:	2.2

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Occlusion spray non flammable green, red, blue**

Revision date: 30.05.2022

Page 11 of 13



Special Provisions: 63, 190, 277, 327, 344, 381, 959  
 Limited quantity: 1000 mL  
 Excepted quantity: E0  
 EmS: F-D, S-U

**Air transport (ICAO-TI/IATA-DGR)**

**14.1. UN number or ID number:** UN 1950  
**14.2. UN proper shipping name:** AEROSOLS, NON-FLAMMABLE  
**14.3. Transport hazard class(es):** 2.2  
**14.4. Packing group:** -  
 Hazard label: 2.2



Special Provisions: A98 A145 A167 A802  
 Limited quantity Passenger: 30 kg G  
 Passenger LQ: Y203  
 Excepted quantity: E0  
 IATA-packing instructions - Passenger: 203  
 IATA-max. quantity - Passenger: 75 kg  
 IATA-packing instructions - Cargo: 203  
 IATA-max. quantity - Cargo: 150 kg

**14.5. Environmental hazards**

ENVIRONMENTALLY HAZARDOUS: No

**14.6. Special precautions for user**

Pressurised container: May burst if heated. Vapours can form explosive mixtures with air.

**14.7. Maritime transport in bulk according to IMO instruments**

not applicable

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulatory information**

Restrictions on use (REACH, annex XVII):

Entry 3, Entry 40, Entry 75

2010/75/EU (VOC): &lt; 95 %

Information according to 2012/18/EU (SEVESO III): Not subject to 2012/18/EU (SEVESO III)

**National regulatory information**

Employment restrictions: Observe restrictions to employment for juveniles according to the 'juvenile work protection guideline' (94/33/EC).

Water hazard class (D): 2 - obviously hazardous to water

**15.2. Chemical safety assessment**

Chemical safety assessments for substances in this mixture were not carried out.

**SECTION 16: Other information****Changes**

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Occlusion spray non flammable green, red, blue**

Revision date: 30.05.2022

Page 12 of 13

This data sheet contains changes from the previous version in section(s): 1,2,4,5,6,7,8,9,10,11,12,13,14,15,16.

**Abbreviations and acronyms**

CLP: Classification, labelling and Packaging  
 REACH: Registration, Evaluation and Authorization of Chemicals  
 GHS: Globally Harmonised System of Classification, Labelling and Packaging of Chemicals  
 UN: United Nations  
 CAS: Chemical Abstracts Service  
 DNEL: Derived No Effect Level  
 DMEL: Derived Minimal Effect Level  
 PNEC: Predicted No Effect Concentration  
 ATE: Acute toxicity estimate  
 LC50: Lethal concentration, 50%  
 LD50: Lethal dose, 50%  
 LL50: Lethal loading, 50%  
 EL50: Effect loading, 50%  
 EC50: Effective Concentration 50%  
 ErC50: Effective Concentration 50%, growth rate  
 NOEC: No Observed Effect Concentration  
 BCF: Bio-concentration factor  
 PBT: persistent, bioaccumulative, toxic  
 vPvB: very persistent, very bioaccumulative  
 ADR: Accord européen sur le transport des marchandises dangereuses par Route  
 (European Agreement concerning the International Carriage of Dangerous Goods by Road)  
 RID: Regulations concerning the international carriage of dangerous goods by rail  
 ADN: European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways  
 (Accord européen relatif au transport international des marchandises dangereuses par voies de navigation intérieures)  
 IMDG: International Maritime Code for Dangerous Goods  
 EmS: Emergency Schedules  
 MFAG: Medical First Aid Guide  
 IATA: International Air Transport Association  
 ICAO: International Civil Aviation Organization  
 MARPOL: International Convention for the Prevention of Marine Pollution from Ships  
 IBC: Intermediate Bulk Container  
 VOC: Volatile Organic Compounds  
 SVHC: Substance of Very High Concern  
 For abbreviations and acronyms, see: ECHA Guidance on information requirements and chemical safety assessment, chapter R.20 (Table of terms and abbreviations).

**Classification for mixtures and used evaluation method according to Regulation (EC) No 1272/2008 [CLP]**

Classification	Classification procedure
Aerosol 3; H229	On basis of test data
STOT SE 3; H336	Bridging principle "Aerosols"
Aquatic Chronic 3; H412	Calculation method

**Relevant H and EUH statements (number and full text)**

H225 Highly flammable liquid and vapour.  
 H229 Pressurised container: May burst if heated.  
 H304 May be fatal if swallowed and enters airways.  
 H319 Causes serious eye irritation.  
 H336 May cause drowsiness or dizziness.  
 H411 Toxic to aquatic life with long lasting effects.  
 H412 Harmful to aquatic life with long lasting effects.  
 EUH066 Repeated exposure may cause skin dryness or cracking.

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Occlusion spray non flammable green, red, blue**

Revision date: 30.05.2022

Page 13 of 13

**Further Information**

The information is based on the present level of our knowledge. It does not, however, give assurance of product properties and establishes no contract legal rights. The receiver of our product is singularly responsible for adhering to existing laws and regulations.

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*