

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Orthocryl LC / Orthocryl LC bite plane paste, transp.**

Revision date: 12.03.2024

Product code: 160-401

Page 1 of 8

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1. Product identifier**

Orthocryl LC / Orthocryl LC bite plane paste, transp.

**Further trade names**

Article No.: 160-401-00, 160-402-00, 160-403-00, 160-404-00, 160-405-00, 160-406-00, 160-410-00

UFI: K7A7-AKEV-FD57-RRU4

**1.2. Relevant identified uses of the substance or mixture and uses advised against****Use of the substance/mixture**

Health services.

**1.3. Details of the supplier of the safety data sheet**

Company name:	DENTAURUM GmbH & Co. KG	
Street:	Turnstr. 31	
Place:	D-75228 Ispringen	
Telephone:	+49 7231 803 0	Telefax: +49 7231 803 295
e-mail:	info@dentaurum.de	
Contact person:	Chemie	
Internet:	www.dentaurum.com	

**1.4. Emergency telephone number:**

+497231803-184 7:00-16:15 (Mo-Do) 7:00-13:15 (Fr)  
+497231803-0 16:15-18:00 (Mo-Do) 13:15-18:00 (Fr)

**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture****Regulation (EC) No 1272/2008**

Hazard categories:  
Skin corrosion/irritation: Skin Irrit. 2  
Serious eye damage/eye irritation: Eye Irrit. 2  
Respiratory or skin sensitisation: Skin Sens. 1  
Hazard Statements:  
Causes skin irritation.  
May cause an allergic skin reaction.  
Causes serious eye irritation.

**2.2. Label elements****Regulation (EC) No 1272/2008****Hazard components for labelling**

1,4-Butanediol dimethacrylate

Signal word: Warning

**Pictograms:****Hazard statements**

H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.

**Precautionary statements**

P261	Avoid breathing dust/fume/gas/mist/vapours/spray.
P272	Contaminated work clothing should not be allowed out of the workplace.

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Orthocryl LC / Orthocryl LC bite plane paste, transp.**

Revision date: 12.03.2024

Product code: 160-401

Page 2 of 8

P280	Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313	If eye irritation persists: Get medical advice/attention.

**2.3. Other hazards**

none

**SECTION 3: Composition/information on ingredients**
**3.2. Mixtures**
**Hazardous components**

CAS No	Chemical name	Quantity
	EC No	
	Index No	
	REACH No	
	Classification (Regulation (EC) No 1272/2008)	
25852-47-5	PEG-400-Dimethacrylate	10-25 %
	Skin Irrit. 2, Eye Irrit. 2; H315 H319	
2082-81-7	1,4-Butanediol dimethacrylate	2,5-10 %
	218-218-1	
	Skin Sens. 1, Aquatic Acute 2; H317 H401	
75980-60-8	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	< 2,5 %
	278-355-8	
	015-203-00-X	
	Repr. 2; H361f	

Full text of H and EUH statements: see section 16.

**SECTION 4: First aid measures**
**4.1. Description of first aid measures**
**After inhalation**

Remove casualty to fresh air and keep warm and at rest. Medical treatment necessary.

**After contact with skin**

After contact with skin, wash immediately with: Water and soap.

**After contact with eyes**

In case of contact with eyes flush immediately with plenty of flowing water for 10 to 15 minutes holding eyelids apart and consult an ophthalmologist.

**After ingestion**

Call a physician immediately.

Rinse mouth immediately and drink 1 glass of water.

**4.2. Most important symptoms and effects, both acute and delayed**

No data available

**4.3. Indication of any immediate medical attention and special treatment needed**

No data available

**SECTION 5: Firefighting measures**
**5.1. Extinguishing media**

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Orthocryl LC / Orthocryl LC bite plane paste, transp.**

Revision date: 12.03.2024

Product code: 160-401

Page 3 of 8

**Suitable extinguishing media**

Foam.  
Extinguishing powder.  
Carbon dioxide (CO<sub>2</sub>).  
Water spray.

**Unsuitable extinguishing media**

Full water jet.

**5.2. Special hazards arising from the substance or mixture**

No data available

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures****General advice**

Remove all sources of ignition.  
Provide adequate ventilation.

**6.2. Environmental precautions**

Do not allow to enter into surface water or drains.  
Do not allow to enter into soil/subsoil.

**6.3. Methods and material for containment and cleaning up****Other information**

Absorb with liquid-binding material (sand, diatomaceous earth, acid- or universal binding agents).  
Treat the recovered material as prescribed in the section on waste disposal.

**6.4. Reference to other sections**

This article doesn't contain hazardous substances or mixtures intended to be released under normal or reasonably foreseeable conditions of use.  
Disposal: see section 13  
Safe handling: see section 7  
Personal protection equipment: see section 8

**SECTION 7: Handling and storage****7.1. Precautions for safe handling****Advice on safe handling**

Keep container tightly closed.

**Advice on protection against fire and explosion**

Keep away from sources of ignition - No smoking.

**Advice on general occupational hygiene**

Remove contaminated, saturated clothing immediately.  
Avoid contact with skin, eyes and clothes.  
Do not breathe gas/fumes/vapour/spray.

**7.2. Conditions for safe storage, including any incompatibilities****Requirements for storage rooms and vessels**

Protect against:

**Hints on joint storage**

Do not store together with:

**7.3. Specific end use(s)**

none

**SECTION 8: Exposure controls/personal protection**

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Orthocryl LC / Orthocryl LC bite plane paste, transp.

Revision date: 12.03.2024

Product code: 160-401

Page 4 of 8

#### 8.1. Control parameters

##### Additional advice on limit values

Does not contain substances above concentration limits fixing an occupational exposure limit.

#### 8.2. Exposure controls

##### Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe dust.

##### Individual protection measures, such as personal protective equipment

##### Hand protection

Recommendation:

Dermatril Single-use gloves. NBR (Nitrile rubber). 0,11 mm

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

##### Skin protection

Protective clothing.

##### Respiratory protection

In case of inadequate ventilation wear respiratory protection.

Respiratory protection necessary at:

## SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

Physical state:	viscous	
Colour:	various, depending on coloration	
Odour:	like:	
Melting point/freezing point:		not determined
Boiling point or initial boiling point and boiling range:		251 °C
Flammability:		not determined
Lower explosion limits:		not determined
Upper explosion limits:		not determined
Flash point:		101 °C
Auto-ignition temperature:		400 °C
Decomposition temperature:		not determined
pH-Value:		not determined
Viscosity / kinematic:		not determined
Water solubility:		partially miscible
Solubility in other solvents		
miscible with most organic solvents		
Partition coefficient n-octanol/water:		not determined
Vapour pressure:		not determined
Vapour pressure:		not determined
Density (at 20 °C):		1,135 g/cm <sup>3</sup>
Bulk density:		not determined
Relative vapour density:		not determined

#### 9.2. Other information

##### Information with regard to physical hazard classes

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Orthocryl LC / Orthocryl LC bite plane paste, transp.**

Revision date: 12.03.2024

Product code: 160-401

Page 5 of 8

**Explosive properties**

none

**Sustaining combustion:**

No data available

**Self-ignition temperature**

Solid:

not determined

Gas:

not determined

**Other safety characteristics****Evaporation rate:**

not determined

**Solvent separation test:**

not determined

**Solvent content:**

not determined

**Solid content:**

not determined

**Sublimation point:**

not determined

**Softening point:**

not determined

**Pour point:**

not determined

**Viscosity / dynamic:**

not determined

**Flow time:**

not determined

**Further Information**

Odour threshold: not determined

**SECTION 10: Stability and reactivity****10.1. Reactivity**

No hazardous reaction when handled and stored according to provisions.

**10.2. Chemical stability**

The product is stable under storage at normal ambient temperatures.

**10.3. Possibility of hazardous reactions**

Reacts with : Peroxides

**10.4. Conditions to avoid**

No data available

**10.5. Incompatible materials**

Reacts with : Peroxides

**10.6. Hazardous decomposition products**

none

**SECTION 11: Toxicological information****11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008****Acute toxicity**

LD50: &gt; 5000 mg/kg (OECD 401; 02.1987)

1,4-Butandioldimethacrylat 10120 mg/kg

**Irritation and corrosivity**

Causes skin irritation.

Causes serious eye irritation.

**Sensitising effects**

May produce an allergic reaction.

**STOT-single exposure**

Based on available data, the classification criteria are not met.

**STOT-repeated exposure**

Based on available data, the classification criteria are not met.

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

### Orthocryl LC / Orthocryl LC bite plane paste, transp.

Revision date: 12.03.2024

Product code: 160-401

Page 6 of 8

#### Aspiration hazard

Based on available data, the classification criteria are not met.

#### 11.2. Information on other hazards

##### Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to humans as no components meets the criteria.

### SECTION 12: Ecological information

#### 12.1. Toxicity

No data available

#### 12.2. Persistence and degradability

No data available

#### 12.3. Bioaccumulative potential

No data available

#### Partition coefficient n-octanol/water

CAS No	Chemical name	Log Pow
2082-81-7	1,4-Butanediol dimethacrylate	3,1

#### 12.4. Mobility in soil

No data available

#### 12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.  
not applicable

#### 12.6. Endocrine disrupting properties

This product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms as no components meets the criteria.

#### 12.7. Other adverse effects

- Do not allow to enter into surface water or drains.
- Do not allow to enter into soil/subsoil.
- Do not allow uncontrolled discharge of product into the environment.

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

##### Disposal recommendations

Send to a hazardous waste incinerator facility under observation of official regulations.

##### List of Wastes Code - used product

070104 WASTES FROM ORGANIC CHEMICAL PROCESSES; wastes from the manufacture, formulation, supply and use (MFSU) of basic organic chemicals; other organic solvents, washing liquids and mother liquors; hazardous waste

##### Contaminated packaging

Dispose of waste according to applicable legislation.

### SECTION 14: Transport information

#### Land transport (ADR/RID)

##### 14.1. UN number or ID number:

No dangerous good in sense of this transport regulation.

##### 14.2. UN proper shipping name:

No dangerous good in sense of this transport regulation.

##### 14.3. Transport hazard class(es):

No dangerous good in sense of this transport regulation.

##### 14.4. Packing group:

No dangerous good in sense of this transport regulation.

#### Inland waterways transport (ADN)

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Orthocryl LC / Orthocryl LC bite plane paste, transp.**

Revision date: 12.03.2024

Product code: 160-401

Page 7 of 8

<b>14.1. UN number or ID number:</b>	No dangerous good in sense of this transport regulation.
<b>14.2. UN proper shipping name:</b>	No dangerous good in sense of this transport regulation.
<b>14.3. Transport hazard class(es):</b>	No dangerous good in sense of this transport regulation.
<b>14.4. Packing group:</b>	No dangerous good in sense of this transport regulation.

**Marine transport (IMDG)**

<b>14.1. UN number or ID number:</b>	No dangerous good in sense of this transport regulation.
<b>14.2. UN proper shipping name:</b>	No dangerous good in sense of this transport regulation.
<b>14.3. Transport hazard class(es):</b>	No dangerous good in sense of this transport regulation.
<b>14.4. Packing group:</b>	No dangerous good in sense of this transport regulation.

**Air transport (ICAO-TI/IATA-DGR)**

<b>14.1. UN number or ID number:</b>	No dangerous good in sense of this transport regulation.
<b>14.2. UN proper shipping name:</b>	No dangerous good in sense of this transport regulation.
<b>14.3. Transport hazard class(es):</b>	No dangerous good in sense of this transport regulation.
<b>14.4. Packing group:</b>	No dangerous good in sense of this transport regulation.

**14.5. Environmental hazards**

ENVIRONMENTALLY HAZARDOUS: No

**14.6. Special precautions for user**

not applicable

**14.7. Maritime transport in bulk according to IMO instruments**

not applicable

**Other applicable information**

No dangerous good in sense of these transport regulations.

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulatory information**

Restrictions on use (REACH, annex XVII):

Entry 3, Entry 75

**National regulatory information**

Water hazard class (D): 2 - obviously hazardous to water

**Additional information**

This mixture contains the following substances of very high concern (SVHC) which are included in the Candidate List according to Article 59 of REACH: none

**15.2. Chemical safety assessment**

For this substance a chemical safety assessment has not been carried out.

**SECTION 16: Other information****Changes**

This data sheet contains changes from the previous version in section(s): 1.

**Relevant H and EUH statements (number and full text)**

H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H361f	Suspected of damaging fertility.
H401	Toxic to aquatic life.

**Further Information**

The above information describes exclusively the safety requirements of the product and is based on our

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**Orthocryl LC / Orthocryl LC bite plane paste, transp.**

Revision date: 12.03.2024

Product code: 160-401

Page 8 of 8

present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*