

**Trade name:** VarseoSmile Crown plus (REF 41107)

**Product no.:** REF 41107

**Current version :** 7.0.0, issued: 17.03.2023

**Replaced version:** 6.1.0, issued: 13.10.2021

**Region:** GB

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1 Product identifier

Trade name

**VarseoSmile Crown plus (REF 41107)**

**VarseoSmile Crown plus A1 (REF 41107 & 41117)**

**VarseoSmile Crown plus A2 (REF 41108 & 41118)**

**VarseoSmile Crown plus A3 (REF 41109 & 41119)**

**VarseoSmile Crown plus B1 (REF 41110 & 41120)**

**VarseoSmile Crown plus B3 (REF 41111 & 41121)**

**VarseoSmile Crown plus BL (REF 41114 & 41124)**

**VarseoSmile Crown plus C2 (REF 41112 & 41122)**

**VarseoSmile Crown plus D3 (REF 41113 & 41123)**

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses of the substance or mixture**

Resin for 3D printing of permanent single crowns, inlays, onlays and veneers

**Uses advised against**

No data available.

### 1.3 Details of the supplier of the safety data sheet

**Address**

BEGO Bremer Goldschlägerei  
Wilh. Herbst GmbH & Co. KG  
Wilhelm-Herbst-Str. 1  
28359 Bremen

Telephone no. +49/ 421/ 2028 – 0  
Fax no. +49/ 421/ 2028 – 115  
e-mail msds@bego.com

**Information provided by / telephone**

Research & Development Department - Materials, alloys and ceramics; +49/ 421/ 2028 – 130 (Chief Development Officer alloys)

**Advice on Safety Data Sheet**

msds@bego.com

### 1.4 Emergency telephone number

For medical advice (in German and English):  
+49 (0)551 192 40 (Giftinformationszentrum Nord)

## SECTION 2: Hazards identification

### 2.1 Classification of the substance or mixture

**Classification in accordance with Regulation (EC) No 1272/2008 (CLP)**

Aquatic Chronic 4; H413

Eye Irrit. 2; H319

Skin Irrit. 2; H315

Skin Sens. 1; H317

STOT SE 3; H335

**Classification information**

**Trade name:** VarseoSmile Crown plus (REF 41107)**Product no.:** REF 41107**Current version :** 7.0.0, issued: 17.03.2023**Replaced version:** 6.1.0, issued: 13.10.2021**Region:** GB

This product is assessed and classified using the methods and criteria below referred to in Article 9 of Regulation (EC) n° 1272/2008:

Physical hazards: determined through assessment data based on the methods or standards referred to in part 2 of Annex I to CLP

Health hazards and environmental hazards: determined through toxicological and ecotoxicological assessment data based on the methods or standards referred to in Part 3, 4 and 5 of Annex I to CLP.

**2.2 Label elements****Labelling according to Regulation (EC) No 1272/2008 (CLP Regulation)****Hazard pictograms**

GHS07

**Signal word**

Warning

**Hazardous component(s) to be indicated on label:**

Esterification products of 4,4'-isopropylidenediphenol, ethoxylated and 2-methylprop-2-enoic acid

**Hazard statement(s)**

H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
H413	May cause long lasting harmful effects to aquatic life.

**Precautionary statement(s)**

P261	Avoid breathing mist/vapours/spray.
P273	Avoid release to the environment.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P312	Call a POISON CENTER/doctor if you feel unwell.
P333+P313	If skin irritation or rash occurs: Get medical advice/attention.
P337+P313	If eye irritation persists: Get medical advice/attention.
P501	Dispose of contents/container to a facility in accordance with local and national regulations.

**2.3 Other hazards**

PBT assessment  
No data available.

vPvB assessment  
No data available.

**SECTION 3: Composition/information on ingredients****3.1 Substances**

Not applicable. The product is not a substance.

**3.2 Mixtures****Hazardous ingredients**

No	Substance name		Additional information	
	CAS / EC / Index / REACH no	Classification (EC) 1272/2008 (CLP)	Concentration	%
1	<b>Esterification products of 4,4'-isopropylidenediphenol, ethoxylated and 2-methylprop-2-enoic acid</b>			
	41637-38-1	Eye Irrit. 2; H319	>= 50.00 - < 75.00	wt%
	-	Skin Irrit. 2; H315		
	-	Skin Sens. 1; H317		
	-	STOT SE 3; H335 Aquatic Chronic 4; H413		
2	<b>diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide</b>		<b>pls. refer to footnote (1)</b>	

**Trade name:** VarseoSmile Crown plus (REF 41107)

**Product no.:** REF 41107

**Current version :** 7.0.0, issued: 17.03.2023

**Replaced version:** 6.1.0, issued: 13.10.2021

**Region:** GB

75980-60-8 278-355-8 015-203-00-X -	Aquatic Chronic 2; H411 Skin Sens. 1B; H317 Repr. 2; H361fd	< 2.50	wt%
--	---	--------	-----

Full Text for all H-phrases and EUH-phrases: pls. see section 16

(1) Aberrant from/in addition to the classification set out in Annex VI, this substance is classified according to European Regulation (EC) No 1272/2008 (CLP), Article 4 (3), paragraph 2.

No	Route, target organ, concrete effect
2	H361fd -; testes, bones; -

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

#### General information

In case of allergic symptoms, especially respiratory tract related, seek medical help immediately. Remove contaminated clothing and shoes immediately, and launder thoroughly before reusing.

#### After inhalation

Remove affected person from the immediate area. Ensure supply of fresh air.

#### After skin contact

In case of contact with skin wash off immediately with copious amounts of water.

#### After eye contact

Remove contact lenses. Rinse eye thoroughly under running water keeping eyelids wide open and protecting the unaffected eye (at least 10 to 15 minutes).

#### After ingestion

Rinse the mouth thoroughly with water. Do not induce vomiting. Never give anything by mouth to an unconscious person.

### 4.2 Most important symptoms and effects, both acute and delayed

No data available.

### 4.3 Indication of any immediate medical attention and special treatment needed

No data available.

## SECTION 5: Firefighting measures

### 5.1 Extinguishing media

#### Suitable extinguishing media

Carbon dioxide; Extinguishing powder; Water spray jet; Alcohol-resistant foam

#### Unsuitable extinguishing media

High power water jet

### 5.2 Special hazards arising from the substance or mixture

In the event of fire, the following can be released: Carbon monoxide (CO); Carbon dioxide (CO<sub>2</sub>); Phosphorus oxides

### 5.3 Advice for firefighters

Use self-contained breathing apparatus. Wear protective clothing.

## SECTION 6: Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

#### For non-emergency personnel

Avoid contact with skin, eyes and clothing. Ensure adequate ventilation. Keep away from ignition sources. Do not inhale vapours. Refer to protective measures listed in sections 7 and 8.

#### For emergency responders

No data available. Personal protective equipment (PPE) - see Section 8.

### 6.2 Environmental precautions

Do not discharge into the drains/surface waters/groundwater. Do not discharge into the subsoil/soil.

**Trade name:** VarseoSmile Crown plus (REF 41107)

**Product no.:** REF 41107

**Current version :** 7.0.0, issued: 17.03.2023

**Replaced version:** 6.1.0, issued: 13.10.2021

**Region:** GB

### 6.3 Methods and material for containment and cleaning up

Take up with absorbent material (e.g., sand, kieselguhr, universal binder). When collected, handle material as described under the section heading "Disposal considerations".

### 6.4 Reference to other sections

No data available.

## SECTION 7: Handling and storage

### 7.1 Precautions for safe handling

#### Advice on safe handling

Provide good ventilation at the work area (local exhaust ventilation, if necessary). Risks inherent to handling the product must be minimised by applying the appropriate protective and preventive measures. Working processes should - so far as possible, according to the state of the art - be designed to rule out bodily contact or the release of hazardous substances.

#### General protective and hygiene measures

Wash hands before breaks and after work. Do not eat, drink or smoke during work time. Remove soiled or soaked clothing immediately. Keep away from food, drink and animal feeding stuffs. Avoid contact with eyes and skin. Provide eye wash fountain in work area.

#### Advice on protection against fire and explosion

Keep away from sources of heat and ignition.

### 7.2 Conditions for safe storage, including any incompatibilities

#### Technical measures and storage conditions

Keep container tightly closed and dry in a cool, well-ventilated place. Protect from heat and direct sunlight.

#### Requirements for storage rooms and vessels

Containers which are opened must be carefully closed and kept upright to prevent leakage. Always keep in containers of same material as the original.

#### Incompatible products

Do not store together with: oxidizing agents

### 7.3 Specific end use(s)

No data available.

## SECTION 8: Exposure controls/personal protection

### 8.1 Control parameters

No parameters available for monitoring.

### 8.2 Exposure controls

#### Appropriate engineering controls

Ensure adequate ventilation, local exhaust at the work station if necessary.

#### Personal protective equipment

##### Respiratory protection

If workplace exposure limits are exceeded, a respiration protection approved for this particular job must be worn. In case of aerosol and mist formation, take appropriate measures for breathing protection in the event workplace threshold values are not specified.

##### Eye / face protection

Safety glasses with side protection shield (EN 166)

##### Hand protection

**Trade name:** VarseoSmile Crown plus (REF 41107)

**Product no.:** REF 41107

**Current version :** 7.0.0, issued: 17.03.2023

**Replaced version:** 6.1.0, issued: 13.10.2021

**Region:** GB

Sufficient protection is given wearing suitable protective gloves checked according to i.e. EN 374, in the event of risk of skin contact with the product. Before use, the protective gloves should be tested in any case for its specific work-station suitability (i.e. mechanical resistance, product compatibility and antistatic properties). Adhere to the manufacturer's instructions and information relating to the use, storage, care and replacement of protective gloves. Protective gloves shall be replaced immediately when physically damaged or worn. Design operations thus to avoid permanent use of protective gloves. Sufficient protection is given wearing suitable protective gloves checked according to i.e. EN 374, in the event of risk of skin contact with the product. Before use, the protective gloves should be tested in any case for its specific work-station suitability (i.e. mechanical resistance, product compatibility and antistatic properties). Adhere to the manufacturer's instructions and information relating to the use, storage, care and replacement of protective gloves. Protective gloves shall be replaced immediately when physically damaged or worn. Design operations thus to avoid permanent use of protective gloves.

Appropriate Material In case of short-term contact / splash protection:

Appropriate Material nitrile  
 Material thickness 0.12 mm  
 Breakthrough time 30 - 60 min

Appropriate Material In case of longer-term contact:

Appropriate Material nitrile  
 Appropriate Material neoprene  
 Material thickness 0.2 mm  
 Breakthrough time 240 - 480 min  
 Appropriate Material nitrile  
 Material thickness >  
 Breakthrough time 480 min

**Other**

Chemical-resistant work clothes.

**Environmental exposure controls**

No data available.

## SECTION 9: Physical and chemical properties

### 9.1 Information on basic physical and chemical properties

<b>State of aggregation</b>	
liquid	
<b>Form</b>	
liquid	
<b>Colour</b>	
Various, depending on coloration	
<b>Odour</b>	
No data available	
<b>pH value</b>	
No data available	
<b>Boiling point / boiling range</b>	
Value	150 °C
<b>Melting point/freezing point</b>	
No data available	
<b>Decomposition temperature</b>	
No data available	
<b>Flash point</b>	
Value	110 °C
<b>Ignition temperature</b>	
No data available	
<b>Flammability</b>	
No data available	

**Trade name:** VarseoSmile Crown plus (REF 41107)

**Product no.:** REF 41107

**Current version :** 7.0.0, issued: 17.03.2023

**Replaced version:** 6.1.0, issued: 13.10.2021

**Region:** GB

Lower explosion limit	
No data available	

Upper explosion limit	
No data available	

Vapour pressure	
No data available	

Relative vapour density	
No data available	

Relative density	
No data available	

Density	
Value	1.4 - 1.5 g/cm <sup>3</sup>
Reference temperature	20 °C

Solubility	
No data available	

Partition coefficient n-octanol/water (log value)			
No	Substance name	CAS no.	EC no.
1	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8
log Pow		3.1	
Reference temperature		23	°C
Method		OECD 117	
Source		ECHA	

Kinematic viscosity	
No data available	

Particle characteristics	
No data available	

## 9.2 Other information

Other information	
No data available.	

## SECTION 10: Stability and reactivity

### 10.1 Reactivity

No data available.

### 10.2 Chemical stability

Stable under recommended storage and handling conditions (See section 7).

### 10.3 Possibility of hazardous reactions

None, when used as directed.

### 10.4 Conditions to avoid

Heat, naked flames and other ignition sources.

### 10.5 Incompatible materials

strong oxidizing agents

### 10.6 Hazardous decomposition products

No hazardous decomposition products known. In case of fire: see section 5.

## SECTION 11: Toxicological information

### 11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute oral toxicity			
No	Substance name	CAS no.	EC no.

Trade name: VarseoSmile Crown plus (REF 41107)

Product no.: REF 41107

Current version : 7.0.0, issued: 17.03.2023

Replaced version: 6.1.0, issued: 13.10.2021

Region: GB

1	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8
LD50	>	5000	mg/kg bodyweight
Species	rat		
Method	OECD 401		
Source	ECHA		

**Acute dermal toxicity**

No	Substance name	CAS no.	EC no.
1	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8
LD50	>	2000	mg/kg bodyweight
Species	rabbit		
Method	OECD 402		
Source	ECHA		

**Acute inhalational toxicity**

No data available

**Skin corrosion/irritation**

No	Substance name	CAS no.	EC no.
1	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8
Species	rabbit		
Source	ECHA		
Evaluation	non-irritant		

**Serious eye damage/irritation**

No	Substance name	CAS no.	EC no.
1	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8
Species	rabbit		
Source	ECHA		
Evaluation	non-irritant		

**Respiratory or skin sensitisation**

No	Substance name	CAS no.	EC no.
1	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8
Route of exposure	Skin		
Species	mouse		
Method	OECD 429		
Source	ECHA		
Evaluation	sensitizing		

**Germ cell mutagenicity**

No	Substance name	CAS no.	EC no.
1	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8
Source	ECHA		
Evaluation/classification	Based on available data, the classification criteria are not met.		

**Reproduction toxicity**

No	Substance name	CAS no.	EC no.
1	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8
Source	ECHA		
Evaluation/classification	Based on available data, the classification criteria are met.		

**Carcinogenicity**

No data available

**STOT - single exposure**

No data available

**STOT - repeated exposure**

No data available

**Aspiration hazard**

No data available

**11.2 Information on other hazards**

Trade name: VarseoSmile Crown plus (REF 41107)

Product no.: REF 41107

Current version : 7.0.0, issued: 17.03.2023

Replaced version: 6.1.0, issued: 13.10.2021

Region: GB

**Endocrine disrupting properties**

No data available.

**Other information**

No data available.

**SECTION 12: Ecological information****12.1 Toxicity**

<b>Toxicity to fish (acute)</b>	
No data available	

<b>Toxicity to fish (chronic)</b>	
No data available	

<b>Toxicity to Daphnia (acute)</b>			
No	Substance name	CAS no.	EC no.
1	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8
EC50		3.53	mg/l
Duration of exposure		48	h
Species	Daphnia magna		
Method	OECD 202		
Source	ECHA		

<b>Toxicity to Daphnia (chronic)</b>	
No data available	

<b>Toxicity to algae (acute)</b>			
No	Substance name	CAS no.	EC no.
1	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8
EC50	>	2.01	mg/l
Duration of exposure		72	h
Species	Pseudokirchneriella subcapitata		
Method	OECD 201		
Source	ECHA		

<b>Toxicity to algae (chronic)</b>	
No data available	

<b>Bacteria toxicity</b>	
No data available	

**12.2 Persistence and degradability**

<b>Biodegradability</b>			
No	Substance name	CAS no.	EC no.
1	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8
Type	aerobic biodegradation		
Value	0	- 10	%
Duration		28	d
Method	OECD 301 F		
Source	ECHA		
Evaluation	not readily biodegradable		

**12.3 Bioaccumulative potential**

<b>Partition coefficient n-octanol/water (log value)</b>			
No	Substance name	CAS no.	EC no.
1	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8
log Pow		3.1	
Reference temperature		23	°C
Method	OECD 117		
Source	ECHA		

**12.4 Mobility in soil**

No data available.



**Trade name:** VarseoSmile Crown plus (REF 41107)

**Product no.:** REF 41107

**Current version :** 7.0.0, issued: 17.03.2023

**Replaced version:** 6.1.0, issued: 13.10.2021

**Region:** GB

## 12.5 Results of PBT and vPvB assessment

Results of PBT and vPvB assessment	
PBT assessment	No data available.
vPvB assessment	No data available.

## 12.6 Endocrine disrupting properties

No data available.

## 12.7 Other adverse effects

No data available.

## 12.8 Other information

Other information
Do not discharge product unmonitored into the environment.

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

#### Product

Allocation of a waste code number, according to the European Waste Catalogue, should be carried out in agreement with the regional waste disposal company.

#### Packaging

Residues must be removed from packaging and when emptied completely disposed of in accordance with the regulations for waste removal. Incompletely emptied packaging must be disposed of in the form of disposal specified by the regional disposer.

## SECTION 14: Transport information

### 14.1 Transport ADR/RID/ADN

The product is not subject to ADR/RID/ADN regulations.

### 14.2 Transport IMDG

The product is not subject to IMDG regulations.

### 14.3 Transport ICAO-TI / IATA

The product is not subject to ICAO-TI / IATA regulations.

### 14.4 Other information

No data available.

### 14.5 Environmental hazards

Information on environmental hazards, if relevant, please see 14.1 - 14.3.

### 14.6 Special precautions for user

No data available.

### 14.7 Maritime transport in bulk according to IMO instruments

Not relevant

## SECTION 15: Regulatory information

### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

#### EU regulations

Regulation (EC) No 1907/2006 (REACH) Annex XIV (List of substances subject to authorisation)
According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances considered as substances requiring authorisation as listed on Annex XIV of the REACH regulation (EC) 1907/2006.

REACH candidate list of substances of very high concern (SVHC) for authorisation

**Trade name:** VarseoSmile Crown plus (REF 41107)**Product no.:** REF 41107**Current version :** 7.0.0, issued: 17.03.2023**Replaced version:** 6.1.0, issued: 13.10.2021**Region:** GB

According to available data and the information provided by preliminary suppliers, the product does not contain substances that are considered substances meeting the criteria for inclusion in annex XIV (List of Substances Subject to Authorisation) as laid down in Article 57 and article 59 of REACH (EC) 1907/2006.

**Regulation (EC) No 1907/2006 (REACH) Annex XVII: RESTRICTIONS ON THE MANUFACTURE, PLACING ON THE MARKET AND USE OF CERTAIN DANGEROUS SUBSTANCES, MIXTURES AND ARTICLES**

The product is considered being subject to REACH regulation (EC) 1907/2006 annex XVII. No 3

The product contains following substance(s) that are considered being subject to REACH regulation (EC) 1907/2006 annex XVII.

No	Substance name	CAS no.	EC no.	No
1	diiron trioxide	1309-37-1	215-168-2	75
2	diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide	75980-60-8	278-355-8	75
3	Iron hydroxide oxide yellow	51274-00-1	257-098-5	75
4	titanium dioxide; [in powder form containing 1 % or more of particles with aerodynamic diameter ≤ 10 µm]	13463-67-7	236-675-5	75
5	triiron-tetraoxide	1317-61-9	215-277-5	75

**Directive 2012/18/EU on the control of major-accident hazards involving dangerous substances**

This product is not subject to Part 1 or 2 of Annex I.

**Other regulations**

Adhere to the national sanitary and occupational safety regulations when using this product.

## 15.2 Chemical safety assessment

A chemical safety assessment has not been carried out for this mixture.

## SECTION 16: Other information

**Sources of key data used to compile the data sheet:**

Regulation (EC) No 1907/2006 (REACH), 1272/2008 (CLP) as amended in each case.

Directives 2000/39/EC, 2006/15/EC, 2009/161/EU, (EU) 2017/164.

National Threshold Limit Values of the corresponding countries as amended in each case.

Transport regulations according to ADR, RID, IMDG, IATA as amended in each case.

The data sources used to determine physical, toxic and ecotoxic data, are indicated directly in the corresponding section.

**Full text of the H- and EUH- phrases drawn up in sections 2 and 3 (provided not already drawn up in these sections)**

H361fd Suspected of damaging fertility. Suspected of damaging the unborn child.

H411 Toxic to aquatic life with long lasting effects.

**Creation of the safety data sheet**

UMCO GmbH - D-21107 Hamburg, Georg-Wilhelm-Strasse 187, Tel.: +49(40)555 546 300, Fax: +49(40)555 546 357, e-mail: umco@umco.de

This information is based on our present knowledge and experience.

The safety data sheet describes products with a view to safety requirements.

It does not however, constitute a guarantee for any specific product properties and shall not establish a legally valid contractual relationship.

**Alterations/supplements:**

Alterations to the previous edition are marked in the left-hand margin.

Document protected by copyright. Alterations or reproductions require the express written permission of UMCO GmbH.

Prod-ID 735423