



Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 1 / GB

Date revised: 09.06.2023

Replaces Version: - / GB

Print date: 09.06.2023

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Drufosoft / Drufosoft pro

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/preparation

Dental splinting material long-term

### 1.3. Details of the supplier of the safety data sheet

#### Address/Manufacturer

Dreve Dentamid GmbH

Max-Planck-Straße 31

59423 Unna

Telephone no. +49 2303 8807-0

Fax no. +49 2303 8807-29

Information provided by / telephone Department Research &amp; Development: Fax: +49 2303 8807-562

E-mail address of

sicherheitsdatenblatt@dreve.com

person responsible

for this SDS

### 1.4. Emergency telephone number

Henkel Fire Department / 24h-Emergency-Contact-No.: +49 211 797-3350

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

This product is not classified hazardous in accordance with Regulation (EC) No 1272/2008.

### 2.2. Label elements

#### Labelling according to regulation (EC) No 1272/2008

#### Supplemental information

EUH210

Safety data sheet available on request.

### 2.3. Other hazards

No special hazards have to be mentioned.

The Substance does not meet PBT-criteria. This substance does not meet the vPvB-criteria. This substance does not have endocrine disrupting properties with respect to humans. This substance does not have endocrine disrupting properties with respect to non-target organisms.

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures



Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 1 / GB

Date revised: 09.06.2023

Replaces Version: - / GB

Print date: 09.06.2023

### Chemical characterization

Ethylene vinyl acetate copolymer (EVA)

### Hazardous ingredients

#### vinyl acetate

CAS No. 108-05-4

EINECS no. 203-545-4

Concentration  $\geq$  0,1 < 1 %

Classification (Regulation (EC) No. 1272/2008)

Flam. Liq. 2 H225

Acute Tox. 4 H332

Carc. 2 H351

STOT SE 3 H335

Additional remarks:

CLP

Regulation (EC) No 1272/2008, Annex VI, Note D

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### General information

No special measures required

#### After inhalation

Ensure supply of fresh air.

#### After skin contact

After contact with molten product cool quickly with cold water. Consult a doctor if skin irritation persists.

#### After eye contact

Separate eyelids, wash the eyes thoroughly with water (15 min.). In case of irritation consult an oculist.

#### After ingestion

Rinse out mouth and give plenty of water to drink.

#### Adhere to personal protective measures when giving first aid

First aider: Pay attention to self-protection!

### 4.2. Most important symptoms and effects, both acute and delayed

Until now no symptoms known so far.

### 4.3. Indication of any immediate medical attention and special treatment needed

#### Hints for the physician / hazards

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

#### Suitable extinguishing media

Recommended: alcohol resistant foam, CO<sub>2</sub>, powders, water spray/mist, Extinguishing measures to suit surroundings

Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 1 / GB

Date revised: 09.06.2023

Replaces Version: - / GB

Print date: 09.06.2023

### 5.3. Advice for firefighters

#### Special protective equipment for fire-fighting

In case of combustion use a suitable breathing apparatus.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

Avoid contact with skin, eyes and clothing.

### 6.2. Environmental precautions

Do not allow to enter drains or waterways.

### 6.3. Methods and material for containment and cleaning up

Pick up mechanically. Clean contaminated floors and objects thoroughly, observing environmental regulations.

### 6.4. Reference to other sections

Refer to protective measures listed in Sections 7 and 8.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

#### Advice on safe handling

Observe the usual precautions for handling chemicals.

#### Advice on protection against fire and explosion

No special measures required.

### 7.2. Conditions for safe storage, including any incompatibilities

#### Requirements for storage rooms and vessels

Store product in closed containers.

#### Hints on storage assembly

Do not store together with foodstuffs.

#### Further information on storage conditions

Keep container tightly closed and dry.

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

#### Other information

Contains no substances with occupational exposure limit values.

### 8.2. Exposure controls

#### General protective and hygiene measures

Observe the usual precautions for handling chemicals.

#### Respiratory protection

Do not inhale vapours

Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 1 / GB

Date revised: 09.06.2023

Replaces Version: - / GB

Print date: 09.06.2023

**Hand protection**

Heat-resistant gloves when using hot material

**Eye protection**

Safety glasses

**Body protection**

Clothing as usual in the chemical industry.

## SECTION 9: Physical and chemical properties

**9.1. Information on basic physical and chemical properties**

<b>Physical state</b>	Thermoforming blank		
<b>Colour</b>	Various, depending on coloration		
<b>Odour</b>	characteristic		
<b>Melting point</b>			
Value	60	to	90 °C
<b>Freezing point</b>			
Remarks	not determined		
<b>Boiling point or initial boiling point and boiling range</b>			
Remarks	not determined		
<b>Flammability</b>			
evaluation	not determined		
<b>Upper and lower explosive limits</b>			
Remarks	not determined		
<b>Flash point</b>			
Remarks	Not applicable		
<b>Ignition temperature</b>			
Remarks	not determined		
<b>Decomposition temperature</b>			
Remarks	not determined		
<b>pH value</b>			
Remarks	not determined		
<b>Viscosity</b>			
Remarks	not determined		
<b>Solubility(ies)</b>			
Remarks	not determined		
<b>Partition coefficient n-octanol/water (log value)</b>			
Remarks	not determined		
<b>Vapour pressure</b>			
Remarks	not determined		
<b>Density and/or relative density</b>			
Value	0,95		g/cm <sup>3</sup>
Temperature	20	°C	
<b>Relative vapour density</b>			
Remarks	not determined		



Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 1 / GB

Date revised: 09.06.2023

Replaces Version: - / GB

Print date: 09.06.2023

**9.2. Other information****Odour threshold**

Remarks not determined

**Evaporation rate (ether = 1) :**

Remarks not determined

**Solubility in water**

Remarks virtually insoluble

**Explosive properties**

evaluation not determined

**Oxidising properties**

Remarks not determined

**Other information**

None known

**SECTION 10: Stability and reactivity****10.1. Reactivity**

No hazardous reactions when stored and handled according to prescribed instructions.

**10.2. Chemical stability**

No hazardous reactions known.

**10.3. Possibility of hazardous reactions**

No hazardous reactions known.

**10.4. Conditions to avoid**

Protect from heat and direct sunlight

**10.5. Incompatible materials**

None known

**SECTION 11: Toxicological information****11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008****Acute oral toxicity**

Remarks Based on available data, the classification criteria are not met.

**Acute oral toxicity (Components)****vinyl acetate**Species rat (male)  
LD50 3470 mg/kg**Acute dermal toxicity**

Remarks Based on available data, the classification criteria are not met.

**Acute dermal toxicity (Components)****vinyl acetate**Species rabbit  
LD50 7440 mg/kg



Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 1 / GB

Date revised: 09.06.2023

Replaces Version: - / GB

Print date: 09.06.2023

**Acute inhalational toxicity**

Remarks Based on available data, the classification criteria are not met.

**Acute inhalative toxicity (Components)****vinyl acetate**

Species	rat		
LC50	4000		ppm(V)
Duration of exposure	4	h	
Administration/Form	Vapors		

**Skin corrosion/irritation**

Remarks Based on available data, the classification criteria are not met.

**Serious eye damage/irritation**

Remarks Based on available data, the classification criteria are not met.

**Sensitization**

Remarks Based on available data, the classification criteria are not met.

**Subacute, subchronic, chronic toxicity**

Remarks not determined

**Mutagenicity**

Remarks Based on available data, the classification criteria are not met.

**Reproductive toxicity**

Remarks Based on available data, the classification criteria are not met.

**Carcinogenicity**

Remarks Based on available data, the classification criteria are not met.

**Carcinogenicity (Components)****vinyl acetate**

evaluation Suspected of causing cancer.

**Specific Target Organ Toxicity (STOT)****Single exposure**

Remarks Based on available data, the classification criteria are not met.

**Repeated exposure**

Remarks Based on available data, the classification criteria are not met.

**Specific Target Organ Toxicity (STOT) (Components)****vinyl acetate****Single exposure**

evaluation May cause respiratory irritation.

**Aspiration hazard**

Based on available data, the classification criteria are not met.

**11.2 Information on other hazards****Endocrine disrupting properties with respect to humans**

This substance does not have endocrine disrupting properties with respect to humans.

**Other information**

No toxicological data are available.

**SECTION 12: Ecological information**

Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 1 / GB

Date revised: 09.06.2023

Replaces Version: - / GB

Print date: 09.06.2023

## 12.1. Toxicity

### General information

not determined

### Fish toxicity (Components)

#### vinyl acetate

Species	Fathead minnow ( <i>Pimephales promelas</i> )		
NOEC	0,551		mg/l
Duration of exposure	34	d	
Method	OECD 210		

### Daphnia toxicity (Components)

#### vinyl acetate

Species	Daphnia magna		
EC50	12,6		mg/l
Duration of exposure	48	h	
Method	OECD 202		

### Algae toxicity (Components)

#### vinyl acetate

Species	Pseudokirchneriella subcapitata		
ErC50	12,7		mg/l
Duration of exposure	72	h	
Method	OECD 201		

## 12.2. Persistence and degradability

### General information

not determined

### Biodegradability (Components)

#### vinyl acetate

Value	82	to	98	%
Duration of test evaluation	14	d		
	Readily biodegradable (according to OECD criteria)			

## 12.3. Bioaccumulative potential

### General information

not determined

### Partition coefficient n-octanol/water (log value)

Remarks not determined

### Octanol/water partition coefficient (log Pow) (Components)

#### vinyl acetate

log Pow	0,73		
Temperature	25	°C	

## 12.4. Mobility in soil

### General information

not determined

## 12.5. Results of PBT and vPvB assessment

### General information

not determined

### Results of PBT and vPvB assessment



Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 1 / GB

Date revised: 09.06.2023

Replaces Version: - / GB

Print date: 09.06.2023

The Substance does not meet PBT-criteria.  
This substance does not meet the vPvB-criteria.

**12.6 Endocrine disrupting properties****Endocrine disrupting properties with respect to the environment**

This substance does not have endocrine disrupting properties with respect to non-target organisms.

**12.7. Other adverse effects****General information**

not determined

**General information / ecology**

Do not discharge product unmonitored into the environment.

**SECTION 13: Disposal considerations****13.1. Waste treatment methods****Disposal recommendations for the product**

Must not be disposed together with household garbage.  
Dispose of waste according to applicable legislation.

**Disposal recommendations for packaging**

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

**SECTION 14: Transport information**

	Land transport ADR/RID	Marine transport IMDG/GGVSee	Air transport ICAO/IATA
<b>14.1. UN number or ID number</b>	The product does not constitute a hazardous substance in land transport.	The product does not constitute a hazardous substance in sea transport.	The product does not constitute a hazardous substance in air transport.
<b>14.2. UN proper shipping name</b>	-	-	-
<b>14.3. Transport hazard class(es)</b>		-	-
Label			
<b>14.4. Packing group</b>		-	-

**SECTION 15: Regulatory information****15.2. Chemical safety assessment**

For this preparation a chemical safety assessment has not been carried out.





Trade name: Drufosoft / Drufosoft pro

Substance number: S0017

Version: 1 / GB

Date revised: 09.06.2023

Replaces Version: - / GB

Print date: 09.06.2023

## SECTION 16: Other information

### Hazard statements listed in Chapter 2/3

H225	Highly flammable liquid and vapour.
H332	Harmful if inhaled.
H335	May cause respiratory irritation.
H351	Suspected of causing cancer.

### CLP categories listed in Chapter 2/3

Acute Tox. 4	Acute toxicity, Category 4
Carc. 2	Carcinogenicity, Category 2
Flam. Liq. 2	Flammable liquid, Category 2
STOT SE 3	Specific target organ toxicity - single exposure, Category 3

### Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: \*\*\*  
This information is based on our present state of knowledge. However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.